



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/14/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:55
In Date: 08/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 07:35
Out Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104574

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4302-01	ETPC5	1	1.12	9.8	9.09	91.8
K4302-02	ETPC6	2	1.14	9.58	9.11	94.4
K4302-03	ETPC7	3	1.15	9.74	9.27	94.5
K4302-04	ETPC8	4	1.14	9.55	9.1	94.6
K4302-05	ETPC9	5	1.14	9.92	8.14	79.7
K4302-06	ETPD0	6	1.14	9.84	8.38	83.2
K4302-07	ETPD1	7	1.15	9.65	8.15	82.4
K4302-08	ETPD2	8	1.14	9.9	8.2	80.6
K4302-09	ETPD3	9	1.13	9.68	7.86	78.7
K4302-10	ETPD4	10	1.12	9.86	8.57	85.2
K4302-11	ETPD5	11	1.13	9.98	8.71	85.6
K4302-12	ETPD6	12	1.12	9.79	7.96	78.9
K4302-13	ETPD7	13	1.12	9.64	8.73	89.3
K4302-14	ETPD8	14	1.12	9.57	8.56	88.0
K4302-15	ETPD9	15	1.12	9.88	8.56	84.9
K4302-16	ETPE0	16	1.13	9.83	8.27	82.1
K4302-17	ETPE1	17	1.13	9.6	8.71	89.5
K4302-18	ETPE2	18	1.15	9.84	8.84	88.5
K4302-19	ETPE2MS	19	1.15	9.84	8.84	88.5
K4302-20	ETPE2MSD	20	1.15	9.84	8.84	88.5
K4302-21	ETPE3	21	1.13	9.97	8.96	88.6

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$