



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/14/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 08:25
In Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: % SOLIDS-OVEN

QC:LB104581

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4326-01	MH1263	1	1.00	2.00	2.00	100.0
K4326-02	MH1264	2	1.00	2.00	2.00	100.0
K4326-03	MH1265	3	1.00	2.00	2.00	100.0
K4326-04	MH1266	4	1.00	2.00	2.00	100.0
K4326-05	MH1267	5	1.00	2.00	2.00	100.0
K4326-06	MH1267D	6	1.00	2.00	2.00	100.0
K4326-07	MH1267S	7	1.00	2.00	2.00	100.0
K4326-08	MH1268	8	1.00	2.00	2.00	100.0
K4326-09	MH1269	9	1.00	2.00	2.00	100.0
K4326-10	MH1270	10	1.00	2.00	2.00	100.0
K4326-11	MH1271	11	1.00	2.00	2.00	100.0
K4326-12	MH1272	12	1.00	2.00	2.00	100.0
K4326-13	MH1273	13	1.00	2.00	2.00	100.0
K4326-14	MH1274	14	1.00	2.00	2.00	100.0
K4326-15	MH1275	15	1.00	2.00	2.00	100.0
K4326-16	MH1276	16	1.00	2.00	2.00	100.0
K4326-17	MH1277	17	1.00	2.00	2.00	100.0
K4326-18	MH1278	18	1.00	2.00	2.00	100.0
K4326-19	MH1279	19	1.00	2.00	2.00	100.0
K4326-20	MH1280	20	1.00	2.00	2.00	100.0
K4326-21	MH1281	21	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$