



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/14/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 08:30
In Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:40
Out Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: % SOLIDS-OVEN

QC:LB104582

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4327-01	MH1282	1	1.00	2.00	2.00	100.0
K4327-02	MH1283	2	1.00	2.00	2.00	100.0
K4327-03	MH1284	3	1.00	2.00	2.00	100.0
K4327-04	MH1285	4	1.00	2.00	2.00	100.0
K4327-05	MH1286	5	1.00	2.00	2.00	100.0
K4327-06	MH1287	6	1.00	2.00	2.00	100.0
K4327-07	MH1288	7	1.00	2.00	2.00	100.0
K4327-08	MH1289	8	1.00	2.00	2.00	100.0
K4327-09	MH1290	9	1.00	2.00	2.00	100.0
K4327-10	MH1291	10	1.00	2.00	2.00	100.0
K4327-11	MH1291D	11	1.00	2.00	2.00	100.0
K4327-12	MH1291S	12	1.00	2.00	2.00	100.0
K4327-13	MH1292	13	1.00	2.00	2.00	100.0
K4327-14	MH1293	14	1.00	2.00	2.00	100.0
K4327-15	MH1294	15	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$