



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 08/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104590

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4264-01	REACTOR-3-C1-C4- (0-2)	1	1.14	9.86	9.02	90.4
K4264-11	REACTOR-3-G1-G4- (0-2)	2	1.14	9.89	8.7	86.4
K4264-16	REACTOR-3-G1-G4- (3)	3	1.13	9.71	8.09	81.1
K4264-21	REACTOR-3-C1-C4- (3)	4	1.14	9.87	8.78	87.5
K4264-22	REACTOR-3-G4- (2)	5	1.14	9.7	8.29	83.5
K4265-01	REACTOR-4-C- (0-3.4)	6	1.15	9.57	8.47	86.9
K4265-04	REACTOR-4-B- (0-3.4)	7	1.14	9.95	8.95	88.6
K4265-12	REACTOR-4-A1-A4- (3.4)	8	1.14	9.55	8.49	87.4
K4265-17	R1-R2- (0-4.3)	9	1.11	9.74	8.73	88.3
K4265-20	REACTOR-4-A1-A4- (0-1)	10	1.14	9.75	8.72	88.0
K4265-21	REACTOR-4-A4- (3.4)	11	1.16	9.91	8.86	88.0
K4267-14	HD-7-4-SE	65	1.12	9.8	7.52	73.7
K4272-01	2S-G1-G4- (0-1)	12	1.12	9.51	8.9	92.7
K4272-06	2S-G1-G4- (1)	13	1.13	9.96	8.95	88.6
K4272-11	2S-K- (0-1)	14	1.15	9.57	7.96	80.9
K4324-01	T8-F1- (0-4)	15	1.15	9.67	8.78	89.6
K4324-02	T8-F1- (0)	16	1.14	9.67	8.72	88.9
K4324-03	T8-F1- (4)	17	1.13	9.88	9.2	92.2
K4324-04	T8-E1- (0-4)	18	1.14	9.84	8.82	88.3
K4324-05	T8-E1- (0)	19	1.13	9.97	9.07	89.8
K4324-06	T8-E1- (4)	20	1.14	9.5	8.4	86.8
K4324-07	T8-C1- (0-4)	21	1.14	9.86	9.09	91.2
K4324-08	T8-C1- (0)	22	1.13	9.73	8.86	89.9
K4324-09	T8-C1- (4)	23	1.13	9.97	9.01	89.1
K4325-02	20190612-LOC-14	24	1.14	9.75	9.19	93.5
K4328-01	20190612-LOC-15-CF	25	1.16	9.9	9.32	93.4
K4329-01	ETGI-292	26	1.14	9.78	7.97	79.1
K4329-02	ETGI-292	27	1.14	9.78	7.97	79.1
K4330-07	0632-COMP-C	28	1.13	9.98	4.7	40.3
K4330-08	0633-COMP-D	29	1.14	9.9	5.42	48.9
K4331-05	0625-COMP-B	30	1.14	9.54	5.4	50.7
K4332-05	0619-COMP-A	31	1.13	9.78	5.15	46.5
K4339-01	T8-A1-A4- (0-4)	32	1.12	9.61	8.56	87.6
K4339-02	T8-A1- (0)	33	1.11	9.64	8.6	87.8
K4339-03	T8-A2- (0)	34	1.12	9.57	8.53	87.7
K4339-04	T8-A3- (0)	35	1.13	9.64	8.71	89.1
K4339-05	T8-A4- (0)	36	1.15	9.65	8.52	86.7
K4339-06	T8-A1- (4)	37	1.13	9.89	8.87	88.4
K4339-07	T8-A2- (4)	38	1.11	9.6	8.58	88.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 08/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104590

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4339-08	T8-A3- (4)	39	1.12	9.5	8.63	89.6
K4339-09	T8-A4- (4)	40	1.13	9.95	8.8	87.0
K4339-10	T8-B1- (0-4)	41	1.11	9.74	8.8	89.1
K4339-11	T8-B1- (0)	42	1.12	9.84	8.93	89.6
K4339-12	T8-B1- (4)	43	1.13	9.96	9.51	94.9
K4340-01	LS-SAMPLE-1	44	1.14	9.96	7.21	68.8
K4340-03	LS-SAMPLE-2	45	1.13	9.6	7.55	75.8
K4341-01	HA-38	46	1.00	2.00	2.00	100.0
K4341-02	HA-39	47	1.00	2.00	2.00	100.0
K4341-03	HA-40	48	1.00	2.00	2.00	100.0
K4341-04	HA-41	49	1.00	2.00	2.00	100.0
K4341-05	HA-42	50	1.00	2.00	2.00	100.0
K4341-06	HA-43	51	1.00	2.00	2.00	100.0
K4341-07	HA-44	52	1.00	2.00	2.00	100.0
K4341-08	HA-45	53	1.00	2.00	2.00	100.0
K4341-09	HA-48	54	1.00	2.00	2.00	100.0
K4342-01	CHRT-26872	55	1.15	9.74	9.32	95.1
K4343-01	OAK-UST-SOIL	56	1.15	9.6	9.51	98.9
K4347-01	HPPER-3-081319-16	57	1.00	2.00	2.00	100.0
K4347-02	HPPER-3-081319-22	58	1.00	2.00	2.00	100.0
K4347-03	HPPER-3-081319-29	59	1.00	2.00	2.00	100.0
K4347-04	HPPER-3-081319-38	60	1.00	2.00	2.00	100.0
K4347-05	HPPER-3-081319-44	61	1.00	2.00	2.00	100.0
K4347-06	HPPER-3-081319-48	62	1.00	2.00	2.00	100.0
K4347-07	HPPER-3-081319-59	63	1.00	2.00	2.00	100.0
K4347-08	HPPER-3-081319-70	64	1.00	2.00	2.00	100.0
K4353-03	40081	66	1.15	9.89	2.74	18.2
K4353-05	30379	67	1.15	9.6	9.17	94.9
K4353-07	40290	68	1.13	9.96	8.59	84.5
K4355-01	RBR-200057-200054	69	1.15	9.55	9.19	95.7
K4356-01	TP-1	70	1.13	9.98	2.14	11.4
K4356-02	TP-1	71	1.13	9.98	2.14	11.4
K4356-03	TP-2	72	1.13	9.54	2.12	11.8
K4356-04	TP-2	73	1.13	9.54	2.12	11.8
K4356-05	TP-3	74	1.15	9.56	2.56	16.8
K4356-06	TP-3	75	1.15	9.56	2.56	16.8
K4356-07	TP-4	76	1.14	9.57	2.23	12.9
K4356-08	TP-4	77	1.14	9.57	2.23	12.9
K4356-09	TP-5	78	1.15	9.56	8.25	84.4



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 08/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104590

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>
K4356-10	TP-5	79	1.15	9.56	8.25	84.4
K4356-11	TP-6	80	1.13	9.88	2.59	16.7
K4356-12	TP-6	81	1.13	9.88	2.59	16.7
K4356-13	TP-7	82	1.14	9.57	2.31	13.9
K4356-14	TP-7	83	1.14	9.57	2.31	13.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$