



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:05
In Date: 08/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104636

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4374-13	MX2H71	1	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-14	MX2H71D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-15	MX2H71S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-16	MX2H72	4	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-17	MX2H72D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-18	MX2H72S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-19	MX2H73	7	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-20	MX2H73D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-21	MX2H73S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-22	MX2H74	10	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-23	MX2H74D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS
K4374-24	MX2H74S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$