



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/19/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:20
Out Date: 08/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104666

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4403-01	T1-H1-H4- (0-2.5)	1	1.12	9.77	8.9	89.9	
K4403-02	T1-H1- (0)	2	1.12	9.52	8.59	88.9	
K4403-03	T1-H2- (0)	3	1.16	9.94	8.93	88.5	
K4403-04	T1-H3- (0)	4	1.16	9.5	8.63	89.6	
K4403-05	T1-H4- (0)	5	1.18	9.85	8.79	87.8	
K4403-06	T1-H1- (2.5)	6	1.15	9.96	9.24	91.8	
K4403-07	T1-H2- (2.5)	7	1.18	9.52	8.56	88.5	
K4403-08	T1-H3- (2.5)	8	1.15	9.97	9.2	91.3	
K4403-09	T1-H4- (2.5)	9	1.18	9.95	8.88	87.8	
K4403-10	T1-G- (0-2.5)	10	1.16	9.77	8.66	87.1	
K4403-11	T1-G- (0)	11	1.18	9.64	8.68	88.7	
K4403-12	T1-G- (2.5)	12	1.14	9.97	8.93	88.2	
K4403-13	T1-I- (0-2.5)	13	1.18	9.88	8.92	89.0	
K4403-14	T1-I- (0)	14	1.17	9.63	8.61	87.9	
K4403-15	T1-I- (2.5)	15	1.18	9.97	8.94	88.3	
K4403-16	T1-J- (0-2.5)	16	1.1	9.8	8.66	86.9	
K4403-17	T1-J- (0)	17	1.16	9.98	8.84	87.1	
K4403-18	T1-J- (2.5)	18	1.13	9.82	8.68	86.9	
K4403-19	T1-K1- (0)	19	1.15	9.64	8.43	85.7	
K4404-01	T1-K1 (2.5)	20	1.17	9.55	8.34	85.6	
K4404-02	T1-L (0-2.5)	21	1.19	9.98	8.94	88.2	
K4404-03	T1-L (0)	22	1.16	9.58	8.46	86.7	
K4404-04	T1-L (2.5)	23	1.14	9.89	8.67	86.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$