



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/20/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 08/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104681

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4095-09	SU-03-081919	1	1.13	9.83	9.35	94.5	
K4095-10	SU-03-081919-1	2	1.13	9.83	9.35	94.5	
K4095-11	SU-03-081919-2	3	1.13	9.83	9.35	94.5	
K4096-04	TR-05-081919	4	1.13	9.81	8.98	90.4	
K4096-05	TR-05-081919	5	1.13	9.81	8.98	90.4	
K4366-01	PX-1	6	1.13	9.81	7.58	74.3	
K4366-02	PX-2	7	1.15	9.98	7.83	75.7	
K4366-03	PX-3	8	1.13	9.58	6.85	67.7	
K4366-04	PX-4	9	1.14	9.93	7.43	71.6	
K4366-05	PX-5	10	1.13	9.52	7.69	78.2	
K4366-06	PX-6	11	1.13	9.98	7.67	73.9	
K4366-07	PX-7	12	1.13	9.9	7.69	74.8	
K4366-08	PX-8	13	1.15	9.82	7.82	76.9	
K4402-01	FRAC-TANK-MUD	14	1.00	2.00	2.00	100.0	
K4406-01	TT-011-IDWSO-20190816	15	1.15	9.71	7.45	73.6	
K4408-01	SP-1	16	1.14	9.67	8.18	82.5	
K4408-02	SP-1A	17	1.13	9.68	8.22	82.9	
K4409-01	S300-S-1.5-2.0-081919	18	1.14	9.77	8.46	84.8	
K4409-10	S301-S-1.5-2.0-081919	19	1.16	9.98	8.74	85.9	
K4409-14	S301-S-1.5-2.0-081919-FD	20	1.14	9.96	8.28	81.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$