



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/21/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:50
Out Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104686

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4411-01	MH1295	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-02	MH1296	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-03	MH1297	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-04	MH1298	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-05	MH1299	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-06	MH12A0	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-07	MH12A1	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-08	MH12A2	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-09	MH12A3	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-10	MH12A4	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-11	MH12A5	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-12	MH12A6	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-13	MH12A7	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-14	MH12A8	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-15	MH12A9	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-16	MH12B0	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-17	MH12B1	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-18	MH12B2	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-19	MH12B3	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-20	MH12B9	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-21	MH12B9D	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4411-22	MH12B9S	22	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$