



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/21/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 08/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104687

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4412-01	MH12B4	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-02	MH12B5	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-03	MH12B6	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-04	MH12B7	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-05	MH12B8	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-06	MH12C0	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-07	MH12C1	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-08	MH12C2	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-09	MH12C3	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-10	MH12C4	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-11	MH12C5	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-12	MH12C6	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-13	MH12C7	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-14	MH12C8	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-15	MH12C8D	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-16	MH12C8S	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-17	MH12C9	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-18	MH12D0	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-19	MH12D1	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-20	MH12D2	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-21	MH12D3	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4412-22	MH12D4	22	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$