



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/21/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:05
In Date: 08/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:41
Out Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104688

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4416-01	ETPZ2	1	1.15	9.7	8.88	90.4	
K4416-02	ETPZ3	2	1.11	9.68	8.95	91.5	
K4416-03	ETPZ4	3	1.12	9.7	8.03	80.5	
K4416-04	ETPZ5	4	1.14	9.81	8.7	87.2	
K4416-05	ETPZ6	5	1.12	9.84	8.58	85.6	
K4416-06	ETPZ6MS	6	1.12	9.84	8.58	85.6	
K4416-07	ETPZ6MSD	7	1.12	9.84	8.58	85.6	
K4416-08	ETPZ7	8	1.13	9.63	8.55	87.3	
K4416-09	ETPZ8	9	1.12	9.89	8.86	88.3	
K4416-10	ETPZ9	10	1.12	9.96	8.64	85.1	
K4416-11	ETQ00	11	1.14	9.92	8.86	87.9	
K4416-12	ETQ01	12	1.11	9.62	8.82	90.6	
K4416-13	ETQ02	13	1.11	9.82	8.87	89.1	
K4416-14	ETQ03	14	1.11	9.97	8.48	83.2	
K4416-15	ETQ04	15	1.13	9.67	8.6	87.5	
K4416-16	ETQ05	16	1.13	9.63	8.81	90.4	
K4416-17	ETQ06	17	1.12	9.6	8.44	86.3	
K4416-18	ETQ07	18	1.11	9.94	8.42	82.8	
K4416-19	ETQ08	19	1.12	9.51	8.17	84.0	
K4416-20	ETQ09	20	1.14	9.93	9.28	92.6	
K4416-21	ETQ10	21	1.13	9.54	8.99	93.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$