



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/22/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104713

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4087-01	HD-02-082119	37	1.12	9.82	8.73	87.5	
K4087-02	HD-02-082119	38	1.12	9.82	8.73	87.5	
K4089-03	OK-01-082119	39	1.15	9.75	9.48	96.9	
K4089-04	OK-01-082119	40	1.13	9.88	9.68	97.7	
K4440-01	S302-J-1.5-2.0-082019	1	1.15	9.83	8.36	83.1	
K4440-02	K4440-01MS	2	1.15	9.83	8.36	83.1	
K4440-03	K4440-01MSD	3	1.15	9.83	8.36	83.1	
K4440-11	S303-J-1.5-2.0-082019	4	1.15	9.96	8.58	84.3	
K4441-02	S59-J-3.0-3.5-082019	5	1.14	9.64	7.98	80.5	
K4441-09	S700-N-0.0-0.5-082019	6	1.19	9.71	8.35	84.0	
K4441-10	K4441-09MS	7	1.19	9.71	8.35	84.0	
K4441-11	K4441-09MSD	8	1.19	9.71	8.35	84.0	
K4443-01	S701-N-0.0-0.5-082019	9	1.16	9.8	7.85	77.4	
K4444-01	SES-GODWIN-PUMP-082019	10	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4444-02	DRIP-PAN-080819	11	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4444-03	DEWATERING-BOX-082019-1	12	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4444-04	DEWATERING-BOX-082019-2	13	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4444-05	DEWATERING-BOX-082019-3	14	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-01	HOPPER-2-082019-16	15	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-02	HOPPER-2-082019-17	16	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-03	HOPPER-2-082019-22	17	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-04	HOPPER-2-082019-29	18	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-05	HOPPER-2-082019-31	19	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-06	HOPPER-2-082019-39	20	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-07	HOPPER-2-082019-56	21	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-08	HOPPER-2-082019-57	22	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-09	HOPPER-2-082019-73	23	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4448-10	HOPPER-2-082019-75	24	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/22/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104713

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4448-11	HOPPER-2-082019-81	25	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-01	MIXING-BOX-1-082019-38	26	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-02	MIXING-BOX-1-082019-46	27	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-03	MIXING-BOX-1-082019-50	28	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-04	MIXING-BOX-1-082019-56	29	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-05	MIXING-BOX-1-082019-58	30	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-06	MIXING-BOX-1-082019-64	31	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-07	MIXING-BOX-1-082019-83	32	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-08	MIXING-BOX-1-082019-97	33	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-09	MIXING-BOX-1-082019-99	34	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-10	MIXING-BOX-1-082019-100	35	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4449-11	MIXING-BOX-1-082019-104	36	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4450-01	WC-A	41	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4450-02	WC-B	42	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4450-03	WC-C	43	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4450-04	WC- (A-C)	44	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4450-05	WC- (A-C) -1	45	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4450-06	WC- (A-C) -2	46	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4450-08	WC-G1	47	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4450-10	WC-G2	48	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4450-12	WC-L1	49	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4450-14	WC-L2	50	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4450-16	WC-L3	51	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4450-18	WC-L4	52	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4450-20	WC-L5	53	1.00	2.00	2.00	100.0	CONCRETE SAMPLE 100 % SOLIDS
K4451-01	B-1	54	1.16	9.64	8.75	89.5	
K4451-02	B-2	55	1.15	9.61	8.89	91.5	
K4451-03	B-3	56	1.14	9.94	9.14	90.9	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/22/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104713

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4452-01	S-CORNER-3	57	1.15	9.94	7.56	72.9	
K4452-02	S-WALL	58	1.15	9.56	7.85	79.7	
K4459-01	HA-09	59	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4459-02	HA-10	60	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4459-03	HA-11	61	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-01	HA-17	62	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-02	HA-18	63	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-03	HA-19	64	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-04	HA-20	65	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-05	HA-21	66	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-06	HA-36	67	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-07	HA-37	68	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-08	HA-38	69	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-09	HA-39	70	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4460-10	HA-40	71	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4461-01	SB-01-9.5-10.0	72	1.14	9.95	9.66	96.7	
K4461-02	SB-01-COMP	73	1.16	9.73	9.39	96.0	
K4461-03	SB-02-9.5-10.0	74	1.18	9.94	9.54	95.4	
K4461-04	SB-02-COMP	75	1.14	9.74	9.28	94.7	
K4461-05	SB-03-10.5-11.0	76	1.14	9.69	9.36	96.1	
K4461-06	SB-03-COMP	77	1.13	9.8	9.52	96.8	
K4461-07	SB-04-11.5-12.0	78	1.13	9.96	8.59	84.5	
K4461-08	SB-04-COMP	79	1.13	9.94	9.56	95.7	
K4461-09	SB-05-8.5-9.0	80	1.13	9.92	9.36	93.6	
K4461-10	SB-05-COMP	81	1.15	9.98	9.8	98.0	
K4461-11	SB-06-7.5-8.0	82	1.12	9.84	9.39	94.8	
K4461-12	SB-06-COMP	83	1.13	9.96	8.76	86.4	
K4461-13	SB-07-8.5-9.0	84	1.12	9.52	9.15	95.6	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/22/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104713

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4461-14	SB-07-COMP	85	1.13	9.86	9.46	95.4	
K4461-15	SB-08-5.5-6.0	86	1.12	9.91	9.56	96.0	
K4461-16	SB-08-COMP	87	1.14	9.85	9.37	94.5	
K4461-17	SB-09-5.5-6.0	88	1.14	9.55	8.7	89.9	
K4461-18	SB-09-COMP	89	1.13	9.69	9.06	92.6	
K4461-19	SB-10-8.5-9.0	90	1.13	9.85	8.29	82.1	
K4461-20	SB-10-COMP	91	1.15	9.64	8.59	87.6	
K4461-21	SB-01-COMP	92	1.16	9.73	9.39	96.0	
K4461-22	SB-02-COMP	93	1.14	9.74	9.28	94.7	
K4461-23	SB-03-COMP	94	1.13	9.8	9.52	96.8	
K4461-24	SB-04-COMP	95	1.13	9.94	9.56	95.7	
K4461-25	SB-05-COMP	96	1.15	9.98	9.8	98.0	
K4461-26	SB-06-COMP	97	1.13	9.96	8.76	86.4	
K4461-27	SB-07-COMP	98	1.13	9.86	9.46	95.4	
K4461-28	SB-08-COMP	99	1.14	9.85	9.37	94.5	
K4461-29	SB-09-COMP	100	1.13	9.69	9.06	92.6	
K4461-30	SB-10-COMP	101	1.15	9.64	8.59	87.6	
K4462-01	SB-46-9.5-10.0	102	1.13	9.95	9.18	91.3	
K4462-02	SB-46-COMP	103	1.13	9.8	9.26	93.8	
K4462-03	SB-47-9.5-10.0	104	1.15	9.7	8.64	87.6	
K4462-04	SB-47-COMP	105	1.11	9.97	9.29	92.3	
K4465-01	FK-PS-01-056	106	1.12	9.83	8.59	85.8	
K4465-02	FK-PS-01-056	107	1.12	9.52	8.47	87.5	
K4468-01	EXT-015-SW085-080119-2	108	1.16	9.98	9.46	94.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$