



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:30
In Date: 08/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:41
Out Date: 08/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104748

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4482-01	ETQ41	1	1.12	9.67	8.62	87.7	
K4482-02	ETQ42	2	1.11	9.55	8.26	84.7	
K4482-03	ETQ43	3	1.12	9.62	8.36	85.2	
K4482-04	ETQ44	4	1.13	9.92	8.83	87.6	
K4482-05	ETQ45	5	1.13	9.75	8.69	87.7	
K4482-06	ETQ46	6	1.12	9.63	8.7	89.1	
K4482-07	ETQ47	7	1.13	9.92	8.8	87.3	
K4482-08	ETQ48	8	1.13	9.7	8.37	84.5	
K4482-09	ETQ49	9	1.14	9.54	8.85	91.8	
K4482-10	ETQ50	10	1.14	9.57	9.00	93.2	
K4482-11	ETQ51	11	1.15	9.59	9.24	95.9	
K4482-12	ETQ52	12	1.14	9.9	9.43	94.6	
K4482-13	ETQ53	13	1.14	9.56	8.72	90.0	
K4482-14	ETQ53MS	14	1.14	9.56	8.72	90.0	
K4482-15	ETQ53MSD	15	1.14	9.56	8.72	90.0	
K4482-16	ETQ54	16	1.15	9.92	8.75	86.7	
K4482-17	ETQ55	17	1.15	9.92	8.7	86.1	
K4482-18	ETQ56	18	1.14	9.61	8.38	85.5	
K4482-19	ETQ57	19	1.14	9.88	8.88	88.6	
K4482-20	ETQ58	20	1.15	9.86	8.64	86.0	
K4482-21	ETQ59	21	1.14	9.71	8.1	81.2	
K4482-22	ETQ60	22	1.14	9.62	8.18	83.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$