



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:14
In Date: 08/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 08/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104769

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4502-01	ETQ75	1	1.15	9.8	8.64	86.6	
K4502-02	ETQ76	2	1.13	9.54	8.61	88.9	
K4502-03	ETQ77	3	1.13	9.93	8.76	86.7	
K4502-04	ETQ78	4	1.14	9.56	8.22	84.1	
K4502-05	ETQ78MS	5	1.14	9.56	8.22	84.1	
K4502-06	ETQ78MSD	6	1.14	9.56	8.22	84.1	
K4502-07	ETQ79	7	1.15	9.85	8.46	84.0	
K4502-08	ETQ80	8	1.15	9.85	8.52	84.7	
K4502-09	ETQ81	9	1.15	9.9	8.58	84.9	
K4502-10	ETQ82	10	1.14	9.8	8.51	85.1	
K4502-11	ETQ83	11	1.14	9.58	8.78	90.5	
K4502-12	ETQ84	12	1.12	9.7	8.76	89.0	
K4502-13	ETQ85	13	1.13	9.72	8.86	90.0	
K4502-14	ETQ86	14	1.15	9.97	8.39	82.1	
K4502-15	ETQ87	15	1.13	9.73	8.66	87.6	
K4502-16	ETQ88	16	1.13	9.66	8.85	90.5	
K4502-17	ETQ89	17	1.12	9.61	9.00	92.8	
K4502-18	ETQ90	18	1.12	9.6	9.01	93.0	
K4502-19	ETQ91	19	1.12	9.63	8.94	91.9	
K4502-20	ETQ92	20	1.12	9.6	8.56	87.7	
K4502-21	ETQ93	21	1.12	9.63	8.45	86.1	
K4502-22	ETQ94	22	1.13	9.6	8.32	84.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$