



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:08
In Date: 08/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 08/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104775

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4409-11	S301-J-2.0-2.5-081919	47	1.15	9.76	8.07	80.4	
K4440-04	S302-J-2.0-2.5-082019	48	1.12	9.92	8.27	81.2	
K4440-05	S302-J-2.0-2.5-082019-F D	49	1.13	9.6	8.01	81.2	
K4441-12	S700-N-1.0-1.5-082019	50	1.12	9.77	8.59	86.4	
K4443-02	S701-N-1.0-1.5-082019	51	1.13	9.81	8.88	89.3	
K4495-01	N1-SB1-1.0-1.5-082219	1	1.15	9.96	8.93	88.3	
K4495-10	N1-SB2-1.0-1.5-082219	2	1.16	9.97	8.76	86.3	
K4495-13	S705-N-0.0-0.5-082219	3	1.12	9.67	8.01	80.6	
K4496-01	S600-M-1.0-1.5-082219	4	1.13	9.93	9.37	93.6	
K4496-04	S600-M-1.0-1.5-082219-F D	5	1.13	9.63	8.87	91.1	
K4496-11	S601-M-1.0-1.5-082219	6	1.12	9.72	8.95	91.0	
K4498-01	SIDEWALL-1	7	1.13	9.93	7.99	78.0	
K4498-02	SIDEWALL-2	8	1.13	9.67	7.72	77.2	
K4498-03	SIDEWALL-3	9	1.15	9.98	7.55	72.5	
K4498-04	BOTTOM-3	10	1.13	9.65	7.67	76.8	
K4498-05	BOTTOM-4	11	1.14	9.76	7.41	72.7	
K4499-01	24408-26706	12	1.14	9.97	8.98	88.8	
K4509-01	SB-48-7.5-8.0	13	1.15	9.65	8.36	84.8	
K4509-02	SB-48-COMP	14	1.14	9.82	8.98	90.3	
K4509-03	SB-49-10.5-11.0	15	1.14	9.66	9.03	92.6	
K4509-04	SB-49-COMP	16	1.14	9.59	9.06	93.7	
K4509-05	SB-50-10.5-11.0	17	1.13	9.65	8.62	87.9	
K4509-06	SB-50-COMP	18	1.13	9.98	9.06	89.6	
K4509-07	SB-51-10.5-11.0	19	1.13	9.97	8.27	80.8	
K4509-08	SB-51-COMP	20	1.16	9.57	8.41	86.2	
K4509-09	SB-45-7.5-8.0	21	1.13	9.96	8.45	82.9	
K4509-10	SB-45-COMP	22	1.14	9.98	8.51	83.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:08
In Date: 08/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 08/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104775

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4510-01	SB-44-6.5-7.0	23	1.13	9.64	6.95	68.4	
K4510-02	SB-44-COMP	24	1.14	9.89	8.92	88.9	
K4510-03	SB-41-11.5-12.0	25	1.14	9.97	8.73	86.0	
K4510-04	SB-41-COMP	26	1.12	9.89	8.26	81.4	
K4510-05	SB-42-9.5-10.0	27	1.14	9.89	8.37	82.6	
K4510-06	SB-42-COMP	28	1.14	9.97	8.41	82.3	
K4510-07	SB-43-8.5-9.0	29	1.14	9.61	8.85	91.0	
K4510-08	SB-43-COMP	30	1.12	9.9	8.65	85.8	
K4510-09	SB-14-7.5-8.0	31	1.14	9.96	8.91	88.1	
K4510-10	SB-14-COMP	32	1.15	9.59	8.45	86.5	
K4510-11	SB-13-6.5-7.0	33	1.15	9.77	8.92	90.1	
K4510-12	SB-13-COMP	34	1.14	9.85	8.47	84.2	
K4510-13	SB-11-10.5-11.0	35	1.14	9.61	9.00	92.8	
K4510-14	SB-11-COMP	36	1.16	9.99	9.06	89.5	
K4513-05	SVOC-GPC-BLANK	37	1.00	2.00	2.00	100.0	
K4513-06	PEST-GPC-BLANK	38	1.00	2.00	2.00	100.0	
K4513-07	PEST-GPC-BLANK-SPIKE	39	1.00	2.00	2.00	100.0	
K4513-08	PCB-GPC-BLANK	40	1.00	2.00	2.00	100.0	
K4513-09	PCB-GPC-BLANK-SPIKE	41	1.00	2.00	2.00	100.0	
K4513-10	SVOC-GPC2-BLANK	42	1.00	2.00	2.00	100.0	
K4513-11	PEST-GPC2-BLANK	43	1.00	2.00	2.00	100.0	
K4513-12	PEST-GPC2-BLANK-SPIKE	44	1.00	2.00	2.00	100.0	
K4513-13	PCB-GPC2-BLANK	45	1.00	2.00	2.00	100.0	
K4513-14	PCB-GPC2-BLANK -SPIKE	46	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$