



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 19:12
In Date: 08/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104794

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4086-05	NB-07-082619	101	1.15	9.92	9.63	96.7	
K4086-06	NB-07-082619-1	102	1.15	9.92	9.63	96.7	
K4086-07	NB-07-082619-2	103	1.15	9.92	9.63	96.7	
K4088-08	CL-02-082619	83	1.13	9.8	9.15	92.5	
K4088-09	CL-02-082619-1	84	1.12	9.68	8.82	90.0	
K4088-10	CL-02-082619-2	85	1.12	9.87	9.36	94.2	
K4088-11	CL-02-082619-3	86	1.15	9.82	9.25	93.4	
K4088-12	CL-02-082619-4	87	1.13	9.56	8.79	90.9	
K4088-13	CL-02-082619-5	88	1.13	9.74	8.96	90.9	
K4092-07	EO-02-082619	104	1.13	9.84	8.92	89.4	
K4092-08	EO-02-082619-1	105	1.13	9.84	8.92	89.4	
K4092-09	EO-02-082619-2	106	1.13	9.84	8.92	89.4	
K4092-10	EO-02-082619-3	107	1.13	9.84	8.92	89.4	
K4093-03	PL-02-082619	108	1.13	9.89	9.41	94.5	
K4093-04	PL-02-082619-1	109	1.13	9.89	9.41	94.5	
K4093-05	PL-02-082619-2	110	1.13	9.89	9.41	94.5	
K4093-06	PL-02-082619-3	111	1.13	9.89	9.41	94.5	
K4093-07	PL-02-082619-4	112	1.13	9.89	9.41	94.5	
K4093-08	PL-02-082619-5	113	1.13	9.89	9.41	94.5	
K4093-09	PL-02-082619-6	114	1.13	9.89	9.41	94.5	
K4506-01	RO-COMP	27	1.16	9.65	4.86	43.6	
K4514-02	S602-M-1.0-1.5-082319	1	1.14	9.82	8.82	88.5	
K4514-03	K4514-02MS	2	1.14	9.82	8.82	88.5	
K4514-04	K4514-02MSD	3	1.14	9.82	8.82	88.5	
K4514-13	S301A-J-1.5-2.0-082319	4	1.14	9.71	8.17	82.0	
K4515-02	S302B-J-1.5-2.0-082319	5	1.14	9.81	8.22	81.7	
K4515-09	S302A-J-1.5-2.0-082319	6	1.13	9.6	7.56	75.9	
K4516-01	T1-S-B2- (0)	7	1.14	9.79	8.72	87.6	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 19:12
In Date: 08/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104794

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4516-02	T1-S-B2- (2.5)	8	1.14	9.71	8.51	86.0	
K4516-03	T1-S-B3- (0)	9	1.14	9.54	8.13	83.2	
K4516-04	T1-S-B3- (2.5)	10	1.15	9.66	8.23	83.2	
K4516-05	T1-S-B4- (0)	11	1.14	9.98	8.63	84.7	
K4516-06	T1-S-B4- (2.5)	12	1.13	9.68	8.18	82.5	
K4516-07	T1-S-B1-B4	13	1.14	9.68	8.42	85.2	
K4516-08	T1-S-E- (0)	14	1.14	9.6	8.22	83.7	
K4516-09	T1-S-E- (2.5)	15	1.13	9.69	8.26	83.3	
K4516-10	T1-S-S-E	16	1.14	9.58	8.23	84.0	
K4516-11	T2-S-B2- (0)	17	1.14	9.66	8.09	81.6	
K4516-12	T2-S-B2- (2.5)	18	1.15	9.7	8.00	80.1	
K4516-13	T2-S-B3- (0)	19	1.14	9.86	8.48	84.2	
K4516-14	T2-S-B3- (2.5)	20	1.14	9.92	8.56	84.5	
K4516-15	T2-S-B4- (0)	21	1.13	9.69	7.92	79.3	
K4516-16	T2-S-B4- (2.5)	22	1.14	9.61	8.15	82.8	
K4516-17	T2-S-B1-B4	23	1.13	9.74	8.01	79.9	
K4516-18	T2-S-D- (0)	24	1.14	9.85	8.44	83.8	
K4516-19	T2-S-D- (2.5)	25	1.15	5.06	4.42	83.6	
K4516-20	T2-S-D	26	1.14	9.58	8.06	82.0	
K4521-01	PX-04	115	1.15	9.85	8.65	86.2	
K4521-02	PX-05	116	1.13	9.65	8.69	88.7	
K4521-03	PX-06	117	1.12	9.58	7.82	79.2	
K4521-04	PX-07	118	1.13	9.77	9.35	95.1	
K4521-05	PX-08	119	1.14	9.9	8.95	89.2	
K4521-06	PX-09	120	1.15	9.85	9.12	91.6	
K4521-07	PX-10	121	1.15	9.97	9.26	92.0	
K4521-08	PX-11	122	1.15	9.67	8.52	86.5	
K4521-09	PX-12	123	1.17	9.56	7.46	75.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 19:12
In Date: 08/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104794

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4521-10	PX-13	124	1.16	9.57	8.92	92.3	
K4521-11	PX-14	125	1.15	9.58	8.88	91.7	
K4521-12	PX-15	126	1.16	9.84	9.22	92.9	
K4521-13	PX-16	127	1.15	9.63	8.07	81.6	
K4521-14	PX-17	128	1.15	9.67	6.52	63.0	
K4521-15	PX-18	129	1.15	9.66	8.42	85.4	
K4524-01	TT-014-IDWSO-20190823	28	1.16	9.94	7.04	67.0	
K4524-02	TT-012-IDWSO-20190823	29	1.15	9.63	8.99	92.5	
K4524-03	TT-013-IDWSO-20190823	30	1.15	9.64	7.79	78.2	
K4525-01	T9-12-13-A1-A4 (0-1.9)	31	1.15	9.76	8.54	85.8	
K4525-02	T9-12-13-A1 (0)	32	1.14	9.78	8.58	86.1	
K4525-03	T9-12-13-A2 (0)	33	1.14	9.44	8.18	84.8	
K4525-04	T9-12-13-A3 (0)	34	1.14	9.75	8.59	86.5	
K4525-05	T9-12-13-A4 (0)	35	1.14	6.92	6.18	87.2	
K4525-06	T9-12-13-A1 (1.9)	36	1.13	9.47	8.28	85.7	
K4525-07	T9-12-13-A2 (1.9)	37	1.15	8.64	8.43	97.2	
K4525-08	T9-12-13-A3 (1.9)	38	1.13	9.28	9.12	98.0	
K4525-09	T9-12-13-A4 (1.9)	39	1.15	9.62	8.59	87.8	
K4525-10	T9-12-13-B1-B4 (0-1.9)	40	1.15	9.95	8.73	86.1	
K4525-11	T9-12-13-B1 (0)	41	1.14	9.9	8.76	87.0	
K4525-12	T9-12-13-B2 (0)	42	1.14	9.94	9.39	93.8	
K4525-13	T9-12-13-B3 (0)	43	1.14	9.66	9.06	93.0	
K4525-14	T9-12-13-B4 (0)	44	1.14	9.93	9.41	94.1	
K4525-15	T9-12-13-B1 (1.9)	45	1.13	9.81	8.41	83.9	
K4525-16	T9-12-13-B2 (1.9)	46	1.13	9.97	8.3	81.1	
K4525-17	T9-12-13-B3 (1.9)	47	1.15	9.56	9.04	93.8	
K4525-18	T9-12-13-B4 (1.9)	48	1.15	9.88	9.3	93.4	
K4525-19	T9-12-13-C (0-1.9)	49	1.13	9.77	8.82	89.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 19:12
In Date: 08/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104794

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4525-20	T9-12-13-C (0)	50	1.14	9.78	9.36	95.1	
K4528-01	T9-12-13-Q2 (0)	51	1.15	9.81	9.18	92.7	
K4528-02	T9-12-13-Q3 (0)	52	1.15	9.9	9.33	93.5	
K4528-03	T9-12-13-Q4 (0)	53	1.16	9.54	8.97	93.2	
K4528-04	T9-12-13-Q2 (1.9)	54	1.15	9.89	9.69	97.7	
K4528-05	T9-12-13-Q3 (1.9)	55	1.16	9.66	8.92	91.3	
K4528-06	T9-12-13-Q4 (1.9)	56	1.16	9.88	9.25	92.8	
K4528-07	T9-12-13-R (0-1.9)	57	1.16	9.73	9.11	92.8	
K4528-08	T9-12-13-R (0)	58	1.16	9.74	9.43	96.4	
K4528-09	T9-12-13-R (1.9)	59	1.16	9.74	9.48	97.0	
K4528-10	T9-12-13-S1-S4 (0-1.9)	60	1.14	9.51	8.56	88.6	
K4528-11	T9-12-13-S1 (0)	61	1.16	9.51	9.18	96.0	
K4528-12	T9-12-13-S2 (0)	62	1.13	9.98	9.03	89.3	
K4528-13	T9-12-13-S3 (0)	63	1.16	9.55	8.58	88.4	
K4528-14	T9-12-13-S4 (0)	64	1.16	9.9	9.61	96.7	
K4528-15	T9-12-13-S1 (1.9)	65	1.14	9.94	9.00	89.3	
K4528-16	T9-12-13-S2 (1.9)	66	1.12	9.98	9.45	94.0	
K4528-17	T9-12-13-S3 (1.9)	67	1.14	9.7	8.93	91.0	
K4528-18	T9-12-13-S4 (1.9)	68	1.12	8.51	7.8	90.4	
K4528-19	T9-12-13-T1-T4 (0-1.9)	69	1.14	9.66	8.88	90.8	
K4529-01	T9-12-13-T2- (0)	130	1.15	9.25	9.06	97.7	
K4529-02	T9-12-13-T3- (0)	131	1.14	9.8	9.02	91.0	
K4529-03	T9-12-13-T4- (0)	132	1.13	8.47	8.22	96.6	
K4529-04	T9-12-13-T1- (1.9)	133	1.13	9.98	8.39	82.0	
K4529-05	T9-12-13-T2- (1.9)	134	1.14	9.98	8.51	83.4	
K4529-06	T9-12-13-T3- (1.9)	135	1.14	9.93	8.79	87.0	
K4529-07	T9-12-13-T4- (1.9)	136	1.14	9.79	9.1	92.0	
K4530-01	390821755.1	70	1.13	9.59	1.42	3.4	SLUDGE SAMPLE



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 19:12
In Date: 08/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104794

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4536-01	SB-1-15-17	71	1.14	9.64	4.37	38.0	
K4536-02	SB-1-15-17	137	1.12	9.51	4.28	37.7	
K4536-03	SB-1-20-22	72	1.16	9.64	4.83	43.3	
K4536-04	SB-1-20-22	138	1.13	7.07	3.41	38.4	
K4536-05	SB-1-25-27	73	1.13	9.94	6.35	59.3	
K4536-06	SB-1-25-27	139	1.13	9.57	6.42	62.7	
K4536-07	SB-1-30-32	74	1.13	9.58	6.42	62.6	
K4536-08	SB-1-30-32	140	1.13	9.92	6.18	57.5	
K4538-01	SB-1	75	1.14	9.8	7.41	72.4	
K4538-02	SB-2	76	1.13	9.75	7.75	76.8	
K4538-03	SB-3	77	1.13	9.85	8.33	82.6	
K4538-04	SB-4	78	1.13	9.88	8.19	80.7	
K4538-05	SB-5	79	1.15	9.6	7.83	79.1	
K4538-06	SB-6	80	1.12	9.71	7.7	76.6	
K4538-07	SB-7	81	1.15	9.8	7.98	79.0	
K4539-01	SB-1	82	1.15	9.69	8.35	84.3	
K4542-01	RUN-228 (IFR) -39	141	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-02	RUN-228 (IFR) -40	142	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-03	RUN-228 (IFR) -41	143	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-04	RUN-228 (IFR) -42	144	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-05	RUN-228 (IFR) -43	145	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-06	RUN-228 (IFR) -44	146	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-07	RUN-228 (IFR) -45	147	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-08	RUN-228 (IFR) -46	148	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-09	RUN-228 (IFR) -47	149	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4542-10	RUN-228 (IFR) -48	150	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4543-01	TP-2-S1	89	1.13	9.77	8.59	86.3	
K4543-02	TP-2-S2	90	1.14	9.77	8.53	85.6	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 19:12
In Date: 08/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104794

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4543-03	TP-2-S	91	1.16	9.88	8.82	87.8	
K4543-05	TP-2-D1	92	1.15	9.52	8.46	87.3	
K4543-06	TP-2-D2	93	1.15	9.57	8.9	92.0	
K4543-07	TP-2-D	94	1.15	9.76	8.38	84.0	
K4543-09	TP-3-S1	95	1.15	9.84	8.82	88.3	
K4543-10	TP-3-S2	96	1.15	9.74	8.35	83.8	
K4543-11	TP-3-S	97	1.15	9.99	9.22	91.3	
K4543-13	TP-9-S1	98	1.13	9.71	9.16	93.6	
K4543-14	TP-9-S2	99	1.14	9.76	9.29	94.5	
K4543-15	TP-9-S	100	1.14	9.74	9.15	93.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$