



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:10
In Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104812

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4371-03	VNJ-222	70	1.13	9.97	9.08	89.9	
K4409-12	S301-S-3.0-3.5-081919	71	1.13	9.77	8.02	79.7	
K4470-02	S702-N-1.0-1.5-082119	72	1.12	9.61	9.02	93.1	
K4470-12	S703-N-1.0-1.5-082119	73	1.14	9.91	8.21	80.6	
K4471-15	S704-N-1.0-1.5-082119	74	1.13	9.77	8.33	83.3	
K4522-02	SIDEWALL-5	75	1.14	9.59	7.25	72.3	
K4526-01	T9-12-13-C (1.9)	1	1.15	9.92	9.51	95.3	
K4526-02	T9-12-13-D (0-1.9)	2	1.14	9.64	9.28	95.8	
K4526-03	T9-12-13-D (0)	3	1.14	9.79	9.42	95.7	
K4526-04	T9-12-13-D (1.9)	4	1.14	9.71	9.34	95.7	
K4526-05	T9-12-13-F1-F4 (0-1.9)	5	1.14	9.94	9.29	92.6	
K4526-06	T9-12-13-F1 (0)	6	1.15	9.62	9.1	93.9	
K4526-07	T9-12-13-F2 (0)	7	1.13	9.77	9.25	94.0	
K4526-08	T9-12-13-F3 (0)	8	1.13	9.65	9.32	96.1	
K4526-09	T9-12-13-F4 (0)	9	1.14	9.53	9.24	96.5	
K4526-10	T9-12-13-F1 (1.9)	10	1.14	9.75	9.06	92.0	
K4526-11	T9-12-13-F2 (1.9)	11	1.14	9.94	9.28	92.5	
K4526-12	T9-12-13-F3 (1.9)	12	1.15	9.68	9.06	92.7	
K4526-13	T9-12-13-F4 (1.9)	13	1.12	9.98	9.53	94.9	
K4526-14	T9-12-13-G1-G4 (0-1.9)	14	1.14	9.82	9.4	95.2	
K4526-15	T9-12-13-G1 (0)	15	1.13	9.97	9.29	92.3	
K4526-16	T9-12-13-G3 (0)	16	1.13	9.7	8.94	91.1	
K4526-17	T9-12-13-G1 (1.9)	17	1.12	9.76	9.04	91.7	
K4526-18	T9-12-13-G3 (1.9)	18	1.12	9.62	9.07	93.5	
K4526-19	T9-12-13-L (0-1.9)	19	1.13	9.82	9.15	92.3	
K4526-20	T9-12-13-L (0)	20	1.13	9.7	9.1	93.0	
K4527-01	T9-12-13-L (1.9)	21	1.12	9.64	8.92	91.5	
K4527-02	T9-12-13-M (0-1.9)	22	1.12	9.95	9.14	90.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:10
In Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104812

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4527-03	T9-12-13-M(0)	23	1.13	9.81	8.93	89.9	
K4527-04	T9-12-13-M(1.9)	24	1.12	9.61	8.83	90.8	
K4527-05	T9-12-13-N1-N4(0-1.9)	25	1.13	9.8	8.78	88.2	
K4527-06	T9-12-13-N1(0)	26	1.12	9.97	8.49	83.3	
K4527-07	T9-12-13-N2(0)	27	1.12	8.96	8.24	90.8	
K4527-08	T9-12-13-N3(0)	28	1.14	9.78	7.76	76.6	
K4527-09	T9-12-13-N4(0)	29	1.13	9.79	9.00	90.9	
K4527-10	T9-12-13-N1(1.9)	30	1.14	9.58	8.09	82.3	
K4527-11	T9-12-13-N2(1.9)	31	1.15	9.91	9.11	90.9	
K4527-12	T9-12-13-N3(1.9)	32	1.14	9.92	8.77	86.9	
K4527-13	T9-12-13-N4(1.9)	33	1.14	9.59	8.84	91.1	
K4527-14	T9-12-13-O(0-1.9)	34	1.14	9.94	8.29	81.3	
K4527-15	T9-12-13-O(0)	35	1.13	9.98	8.19	79.8	
K4527-16	T9-12-13-O(1.9)	36	1.12	9.5	8.26	85.2	
K4527-17	T9-12-13-P(0-1.9)	37	1.14	9.58	9.42	98.1	
K4527-18	T9-12-13-P(0)	38	1.13	9.87	9.56	96.5	
K4527-19	T9-12-13-P(1.9)	39	1.15	9.73	9.54	97.8	
K4527-20	T9-12-13-Q2-Q4(0-1.9)	40	1.14	9.83	9.03	90.8	
K4528-20	T9-12-13-T1(0)	125	1.16	7.57	7.43	97.8	
K4544-01	T16-17-B1-B4-(0-2.25)	83	1.14	9.62	8.39	85.5	
K4544-02	T16-17-B1-(0)	84	1.13	9.7	9.35	95.9	
K4544-03	T16-17-B2-(0)	85	1.12	9.75	9.51	97.2	
K4544-04	T16-17-B3-(0)	86	1.14	9.59	8.86	91.4	
K4544-05	T16-17-B4-(0)	87	1.12	9.67	9.1	93.3	
K4544-06	T16-17-B1-(2.25)	88	1.15	9.81	9.39	95.2	
K4544-07	T16-17-B2-(2.25)	89	1.14	9.8	9.34	94.7	
K4544-08	T16-17-B3-(2.25)	90	1.13	9.83	8.15	80.7	
K4544-09	T16-17-B4-(2.25)	91	1.13	9.71	9.26	94.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:10
In Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104812

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4544-10	T16-17-C1-C4- (0-2.25)	92	1.13	9.81	9.3	94.1	
K4544-11	T16-17-C1- (0)	93	1.13	9.96	9.69	96.9	
K4544-12	T16-17-C2- (0)	94	1.13	9.64	9.48	98.1	
K4544-13	T16-17-C3- (0)	95	1.13	9.87	9.43	95.0	
K4544-14	T16-17-C4- (0)	96	1.15	9.68	9.25	95.0	
K4544-15	T16-17-C1- (2.25)	97	1.15	9.9	9.54	95.9	
K4544-16	T16-17-C2- (2.25)	98	1.14	9.73	9.38	95.9	
K4544-17	T16-17-C3- (2.25)	99	1.15	9.82	9.41	95.3	
K4544-18	T16-17-C4- (2.25)	100	1.16	9.85	9.48	95.7	
K4544-19	T16-17-E1-E4- (0-2.25)	101	1.15	9.58	8.93	92.3	
K4544-20	T16-17-E1- (0)	102	1.13	9.6	9.17	94.9	
K4545-01	T16-17-E2- (0)	126	1.13	9.81	9.08	91.6	
K4545-02	T16-17-E3- (0)	127	1.14	9.5	8.7	90.4	
K4545-03	T16-17-E4- (0)	128	1.16	9.47	8.13	83.9	
K4545-04	T16-17-E1- (2.25)	129	1.16	9.64	9.17	94.5	
K4545-05	T16-17-E2- (2.25)	130	1.16	9.62	9.09	93.7	
K4545-06	T16-17-E3- (2.25)	131	1.15	9.93	9.06	90.1	
K4545-07	T16-17-E4- (2.25)	132	1.13	9.67	9.11	93.4	
K4546-01	T3-S-C1 (0)	41	1.13	9.89	8.77	87.2	
K4546-02	T3-S-C1 (2.5)	42	1.13	9.91	8.46	83.5	
K4546-03	T3-S-C2 (0)	43	1.16	9.65	8.44	85.7	
K4546-04	T3-S-C2 (2.5)	44	1.14	9.83	8.6	85.8	
K4546-05	T3-S-C3 (0)	45	1.14	9.87	8.72	86.8	
K4546-06	T3-S-C3 (2.5)	46	1.15	9.98	8.43	82.4	
K4546-07	T3-S-C4 (0)	47	1.13	9.69	8.31	83.9	
K4546-08	T3-S-C4 (2.5)	48	1.13	9.58	8.13	82.8	
K4546-09	T3-S-C1-C4	49	1.14	9.74	8.29	83.1	
K4546-10	T3-S-A (0)	50	1.13	9.76	8.53	85.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:10
In Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104812

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4546-11	T3-S-A(2.5)	51	1.14	9.58	8.14	82.9	
K4546-12	T3-S-A	52	1.12	9.67	8.27	83.6	
K4546-13	T3-S-B1(0)	53	1.13	9.82	8.36	83.2	
K4546-14	T3-S-B1(2.5)	54	1.13	9.7	8.17	82.1	
K4546-15	T3-S-B2(0)	55	1.12	9.77	8.94	90.4	
K4546-16	T3-S-B2(2.5)	56	1.12	9.96	9.44	94.1	
K4546-17	T3-S-B3(0)	57	1.13	9.89	8.89	88.6	
K4546-18	T3-S-B3(2.5)	58	1.13	9.52	8.58	88.8	
K4546-19	T3-S-B4(0)	59	1.12	9.81	9.16	92.5	
K4546-20	T3-S-B4(2.5)	60	1.12	9.86	9.11	91.4	
K4547-01	T3-S-B1-B4	103	1.15	9.72	8.59	86.8	
K4547-02	T3-S-D(0)	104	1.15	9.81	8.53	85.2	
K4547-03	T3-S-D(2.5)	105	1.15	9.57	9.03	93.6	
K4547-04	T3-S-D	106	1.16	9.91	9.09	90.6	
K4547-05	T3-S-G(0)	107	1.12	9.56	8.65	89.2	
K4547-06	T3-S-G(2.5)	108	1.14	9.67	8.6	87.5	
K4547-07	T3-S-G	109	1.12	9.66	8.78	89.7	
K4547-08	T3-S-F(0)	110	1.14	9.57	8.35	85.5	
K4547-09	T3-S-F(2.5)	111	1.12	9.51	8.64	89.6	
K4547-10	T3-S-F	112	1.14	9.65	8.57	87.3	
K4547-11	T7-S-G1(0)	113	1.12	9.57	8.81	91.0	
K4547-12	T7-S-G1(2.5)	114	1.12	9.63	8.33	84.7	
K4547-13	T7-S-G2(0)	115	1.15	9.72	8.48	85.5	
K4547-14	T7-S-G2(2.5)	116	1.13	9.65	8.47	86.2	
K4547-15	T7-S-G1-G4	117	1.14	9.9	9.65	97.1	
K4547-16	T12-1(0)	118	1.12	9.63	8.26	83.9	
K4547-17	T12-1(2.5)	119	1.14	9.66	8.19	82.7	
K4547-18	T12-2(0)	120	1.13	9.68	8.61	87.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:10
In Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104812

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4547-19	T11-1 (0)	121	1.13	9.67	8.8	89.8	
K4547-20	T11-1 (2.5)	122	1.13	9.93	8.73	86.4	
K4548-01	T11-2 (0)	61	1.13	9.92	8.41	82.8	
K4548-02	T11-2 (2.5)	62	1.14	9.72	7.61	75.4	
K4548-03	T8-S-G1 (0)	63	1.13	9.62	8.37	85.3	
K4548-04	T8-S-G1 (2.5)	64	1.13	9.77	8.54	85.8	
K4548-05	T8-S-G2 (0)	65	1.13	9.89	8.72	86.6	
K4548-06	T8-S-G2 (2.5)	66	1.16	9.97	8.62	84.7	
K4548-07	T8-S-G3 (0)	67	1.13	9.98	8.74	86.0	
K4548-08	T8-S-G3 (2.5)	68	1.12	9.75	8.59	86.6	
K4548-09	T8-S-G1-G4 (0)	69	1.14	9.91	8.57	84.7	
K4549-01	S700A-N-0.0-0.5-082619	76	1.13	9.88	8.25	81.4	
K4551-01	S401-K1-1.0-1.5-082619-FD	77	1.15	9.64	7.32	72.7	
K4551-08	S400-K1-1.0-1.5-082619	78	1.14	9.86	7.76	75.9	
K4551-17	S402-K1-1.0-1.5-082619	79	1.13	9.83	8.7	87.0	
K4551-18	K4551-17MS	80	1.15	9.53	8.63	89.3	
K4551-20	K4551-17MSD	81	1.15	9.81	8.82	88.6	
K4559-01	SOIL-SAMPLE-4	82	1.14	9.76	7.56	74.5	
K4560-01	S404-K1-1.0-1.5-082619	123	1.13	9.95	8.06	78.6	
K4560-10	S401-K1-1.0-1.5-082619	124	1.14	9.84	7.22	69.9	
K4562-01	SB-02-15-17	133	1.13	9.72	5.5	50.9	
K4562-02	SB-02-15-17	134	1.12	9.71	5.29	48.5	
K4562-03	SB-02-20-22	135	1.12	9.8	4.55	39.5	
K4562-04	SB-02-20-22	136	1.16	9.6	4.44	38.9	
K4562-05	SB-02-25-27	137	1.13	9.81	6.48	61.6	
K4562-06	SB-02-25-27	138	1.14	9.89	6.33	59.3	
K4562-07	SB-02-30-32	139	1.13	9.95	6.65	62.6	
K4562-08	SB-02-30-32	140	1.14	9.68	6.47	62.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 18:10
In Date: 08/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104812

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4563-01	CIW-6S- (0-4)	141	1.13	9.72	8.32	83.7	
K4567-01	TP-7-S1	142	1.12	9.59	8.94	92.3	
K4567-02	TP-7-S2	143	1.12	9.6	8.88	91.5	
K4567-03	TP-7-S3	144	1.12	9.65	8.94	91.7	
K4567-04	TP-7-S4	145	1.12	9.51	8.89	92.6	
K4567-05	TP-7-S	146	1.12	9.71	9.08	92.7	
K4567-07	TP-7-D1	147	1.15	9.8	8.86	89.1	
K4567-08	TP-7-D2	148	1.14	9.91	9.17	91.6	
K4567-09	TP-7-D3	149	1.13	9.7	8.8	89.5	
K4567-10	TP-7-D4	150	1.12	9.81	9.13	92.2	
K4567-11	TP-7-D	151	1.12	9.72	8.84	89.8	
K4567-13	TP-8-S1	152	1.14	9.51	8.61	89.2	
K4567-14	TP-8-S2	153	1.13	9.72	8.49	85.7	
K4567-15	TP-8-S	154	1.12	9.79	8.8	88.6	
K4567-17	TP-8-D1	155	1.13	9.66	8.84	90.4	
K4567-18	TP-8-D2	156	1.12	9.98	8.82	86.9	
K4567-19	TP-8-D	157	1.12	9.69	8.91	90.9	
K4571-01	RT-3425	158	1.15	9.76	8.53	85.7	
K4571-02	RT-3425	159	1.15	9.76	8.53	85.7	
K4571-04	RBR-200028	160	1.12	9.84	7.28	70.6	
K4571-05	RBR-200028	161	1.12	9.84	7.28	70.6	
K4571-07	RBR-200052	162	1.12	9.52	7.99	81.8	
K4571-08	RBR-200052	163	1.12	9.52	7.99	81.8	
K4572-01	RBR-200059-200035	164	1.14	9.71	8.6	87.0	
K4576-01	RT-1868-4221	165	1.14	9.75	8.33	83.5	
K4576-02	RBR-200034-RBR-200036	166	1.12	9.88	9.06	90.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$