



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/29/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:20
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104834

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4517-01	T2-S-A1-(0)	1	1.13	9.65	8.12	82.0	
K4517-02	T2-S-A1-(2.5)	2	1.13	9.92	8.67	85.8	
K4517-03	T2-S-A2-(0)	3	1.12	9.71	8.38	84.5	
K4517-04	T2-S-A2-(2.5)	4	1.13	9.9	8.54	84.5	
K4517-05	T2-S-A1-A4	5	1.14	9.71	8.2	82.4	
K4517-08	T4-S-D2-(0)	6	1.15	9.68	8.36	84.5	
K4517-09	T4-S-D2-(2.5)	7	1.13	9.85	9.19	92.4	
K4517-10	T4-S-D4-(0)	8	1.13	9.74	8.31	83.4	
K4517-11	T4-S-D4-(2.5)	9	1.15	9.97	8.61	84.6	
K4517-12	T4-S-D1-D4	10	1.13	9.76	8.74	88.2	
K4517-13	T4-S-C-(0)	11	1.12	9.74	8.8	89.1	
K4517-14	T4-S-C-(2.5)	12	1.13	9.67	8.22	83.0	
K4517-15	T4-S-C	13	1.14	9.93	8.54	84.2	
K4517-16	T4-S-E-(0)	14	1.15	9.61	8.22	83.6	
K4517-17	T4-S-E-(2.5)	15	1.15	9.97	8.49	83.2	
K4517-18	T4-S-E	16	1.13	9.98	8.37	81.8	
K4568-01	SB-12-8.5-9.0	28	1.15	9.81	9.22	93.2	
K4568-02	SB-12-COMP	29	1.14	9.75	9.44	96.4	
K4568-03	SB-15-7.5-8.0	30	1.13	9.64	9.06	93.2	
K4568-04	SB-15-COMP	31	1.14	9.7	9.2	94.2	
K4568-05	SB-16-9.5-10.0	32	1.13	9.78	9.11	92.3	
K4568-06	SB-16-COMP	33	1.14	9.75	8.78	88.7	
K4568-07	SB-17-8.5-9.0	34	1.12	9.63	9.04	93.1	
K4568-08	SB-17-COMP	35	1.13	9.65	9.31	96.0	
K4568-09	SB-18-6.5-7.0	36	1.12	9.88	9.22	92.5	
K4568-10	SB-18-COMP	37	1.12	9.76	9.13	92.7	
K4568-11	SB-19-10-10.5	38	1.12	9.87	8.66	86.2	
K4568-12	SB-19-COMP	39	1.13	9.77	8.57	86.1	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/29/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:20
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104834

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4568-13	SB-37-10.5-11.0	40	1.12	9.65	9.16	94.3	
K4568-14	SB-37-COMP	41	1.13	9.77	9.11	92.4	
K4568-15	SB-36-10.5-11.0	42	1.16	9.69	9.25	94.8	
K4568-16	SB-36-COMP	43	1.13	9.64	9.05	93.1	
K4568-17	SB-35-9.5-10.0	44	1.13	9.6	8.69	89.3	
K4568-18	SB-35-COMP	45	1.14	9.75	8.72	88.0	
K4568-19	SB-38-7.5-8.0	46	1.13	9.66	9.31	95.9	
K4568-20	SB-38-COMP	47	1.13	9.75	9.4	95.9	
K4569-01	SB-39-6.5-7.0	17	1.13	9.82	9.4	95.2	
K4569-02	SB-39-COMP	18	1.12	9.62	8.22	83.5	
K4573-01	S700B-N-0.0-0.5-082719	19	1.14	9.84	9.35	94.4	
K4573-10	S702A-N-0.0-0.5-082719	20	1.13	9.53	9.31	97.4	
K4573-11	K4573-10MS	21	1.13	9.53	9.31	97.4	
K4573-12	K4573-10MSD	22	1.13	9.53	9.31	97.4	
K4574-04	S701A-N-0.0-0.5-082719	23	1.13	9.7	9.11	93.1	
K4574-12	S702B-N-0.0-0.5-082719	24	1.12	9.85	8.49	84.4	
K4574-13	S702B-N-1.0-1.5-082719	25	1.13	9.56	9.35	97.5	
K4575-04	S701B-N-0.0-0.5-082719	26	1.12	9.63	8.22	83.4	
K4575-13	S703A-N-0.0-0.5-082719	27	1.14	9.6	8.19	83.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$