



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/30/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:30
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104842

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4585-01	MH12H9	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-02	MH12J0	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-03	MH12J1	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-04	MH12J2	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-05	MH12J3	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-06	MH12J3D	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-07	MH12J3S	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-08	MH12J5	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-09	MH12J6	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-10	MH12J7	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-11	MH12J8	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-12	MH12J9	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-13	MH12K0	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-14	MH12K1	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-15	MH12K2	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-16	MH12K3	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-17	MH12K4	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-18	MH12K5	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-19	MH12K6	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4585-20	MH12K7	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$