



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/30/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:45
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:38
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104843

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4586-01	MH12J4	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-02	MH12J4D	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-03	MH12J4S	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-04	MH12K8	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-05	MH12K9	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-06	MH12L0	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-07	MH12L1	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-08	MH12L2	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-09	MH12L3	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-10	MH12L4	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-11	MH12L5	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-12	MH12L6	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-13	MH12L7	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-14	MH12L8	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-15	MH12L9	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-16	MH12M0	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-17	MH12M1	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-18	MH12M2	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-19	MH12M3	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4586-20	MH12M4	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$