



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104844

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4095-13	SU-01-082819	79	1.12	9.96	9.4	93.7	
K4095-14	SU-01-082819-1	80	1.13	9.54	9.08	94.5	
K4095-15	SU-01-082819-2	81	1.15	9.6	8.87	91.4	
K4439-02	S-1-(1-3)	58	1.13	9.8	8.68	87.1	
K4439-03	S-2-(1-3)	59	1.13	9.85	9.1	91.4	
K4439-04	S-1-(5-7)	60	1.13	9.68	8.8	89.7	
K4439-05	S-4A-(1-3)	61	1.14	9.7	8.82	89.7	
K4439-06	S-3-(1-3)	62	1.14	9.81	8.79	88.2	
K4439-07	S-5-(1-3)	63	1.12	9.79	8.09	80.4	
K4439-08	S-2-(5-7)	64	1.12	9.76	9.00	91.2	
K4439-09	S-4A-(5-7)	65	1.11	9.76	8.77	88.6	
K4439-10	S-5-(5-7)	66	1.12	9.55	9.03	93.8	
K4441-13	S700-N-2.0-2.5-082019	1	1.14	9.74	8.59	86.6	
K4579-01	PAR-12-A1-A4-(0-3.5)	2	1.13	9.89	8.71	86.5	
K4579-02	PAR-12-A1-(0)	3	1.11	9.93	8.65	85.5	
K4579-03	PAR-12-A3-(0)	4	1.12	9.97	8.9	87.9	
K4579-04	PAR-12-A4-(0)	5	1.12	9.76	8.2	81.9	
K4579-05	PAR-12-A1-(3.5)	6	1.12	9.9	8.77	87.1	
K4579-06	PAR-12-A3-(3.5)	7	1.14	9.75	8.56	86.2	
K4579-07	PAR-12-A4-(3.5)	8	1.12	9.7	8.48	85.8	
K4579-08	PAR-12-B1-B4-(0-3.5)	9	1.12	9.74	8.92	90.5	
K4579-09	PAR-12-B1-(0)	10	1.12	8.5	7.88	91.6	
K4579-10	PAR-12-B2-(0)	11	1.13	9.89	8.98	89.6	
K4579-11	PAR-12-B3-(0)	12	1.16	9.93	8.81	87.2	
K4579-12	PAR-12-B4-(0)	13	1.13	9.85	9.12	91.6	
K4579-13	PAR-12-B1-(3.5)	14	1.14	8.69	8.08	91.9	
K4579-14	PAR-12-B2-(3.5)	15	1.13	9.81	8.98	90.4	
K4579-15	PAR-12-B3-(3.5)	16	1.13	9.82	8.02	79.3	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104844

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4579-16	PAR-12-B4- (3.5)	17	1.12	9.62	8.93	91.9	
K4579-17	PAR-12-C1- (0)	18	1.12	9.81	7.79	76.8	
K4580-01	T12-A3- (3.5)	19	1.12	9.51	9.24	96.8	
K4580-02	T12-A4- (3.5)	20	1.12	9.98	9.54	95.0	
K4580-03	T12-B1-B4- (0-3.5)	21	1.13	9.95	9.02	89.5	
K4580-04	T12-B1- (0)	22	1.14	9.74	9.36	95.6	
K4580-05	T12-B2- (0)	23	1.12	9.79	9.3	94.3	
K4580-06	T12-B3- (0)	24	1.14	9.69	9.01	92.0	
K4580-07	T12-B4- (0)	25	1.13	9.61	7.89	79.7	
K4580-08	T12-B1- (3.5)	26	1.12	9.67	8.74	89.1	
K4580-09	T12-B2- (3.5)	27	1.12	9.87	9.24	92.8	
K4580-10	T12-B3- (3.5)	28	1.13	9.08	8.36	90.9	
K4580-11	PAR-12-C1- (3.5)	29	1.13	9.63	7.73	77.6	
K4580-12	PAR-12-C2- (0)	30	1.12	9.83	8.59	85.8	
K4580-13	PAR-12-C2- (3.5)	31	1.13	9.78	8.77	88.3	
K4580-14	T12-A1-A4- (0-3.5)	32	1.14	9.77	9.02	91.3	
K4580-15	T12-A1- (0)	33	1.12	9.88	9.44	95.0	
K4580-16	T12-A2- (0)	34	1.13	9.5	9.35	98.2	
K4580-17	T12-A3- (0)	35	1.13	9.98	9.5	94.6	
K4580-18	T12-A4- (0)	36	1.13	9.84	9.14	92.0	
K4580-19	T12-A1- (3.5)	37	1.13	9.51	8.73	90.7	
K4580-20	T12-A2- (3.5)	38	1.14	9.81	9.26	93.7	
K4581-01	T12-B4- (3.5)	39	1.11	9.93	8.57	84.6	
K4581-02	T12-C1-C4- (0-3.5)	40	1.12	9.9	9.18	91.8	
K4581-03	T12-C1- (0)	41	1.12	9.72	9.38	96.0	
K4581-04	T12-C2- (0)	42	1.12	9.72	8.5	85.8	
K4581-05	T12-C3- (0)	43	1.13	9.52	8.7	90.2	
K4581-06	T12-C4- (0)	44	1.15	9.67	8.81	89.9	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 08/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104844

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4581-07	T12-C1-(3.5)	45	1.14	9.79	9.04	91.3	
K4581-08	T12-C2-(3.5)	46	1.13	9.65	8.27	83.8	
K4581-09	T12-C3-(3.5)	47	1.14	9.98	9.54	95.0	
K4581-10	T12-C4-(3.5)	48	1.15	9.83	8.36	83.1	
K4581-11	T12-D1-D4-(0-3.5)	49	1.13	9.58	8.95	92.5	
K4581-12	T12-D1-(0)	50	1.12	9.75	9.04	91.8	
K4581-13	T12-D2-(0)	51	1.12	9.98	8.75	86.1	
K4581-14	T12-D3-(0)	52	1.14	9.97	9.61	95.9	
K4581-15	T12-D4-(0)	53	1.12	9.84	9.47	95.8	
K4581-16	T12-D1-(3.5)	54	1.14	9.8	8.75	87.9	
K4581-17	T12-D2-(3.5)	55	1.12	9.61	8.83	90.8	
K4581-18	T12-D3-(3.5)	56	1.14	9.98	8.89	87.7	
K4581-19	T12-D4-(3.5)	57	1.14	9.6	9.03	93.3	
K4582-01	RT-3422	70	1.14	9.62	8.62	88.2	
K4583-01	GS-STOCKPILE	71	1.14	9.55	9.08	94.4	
K4587-01	FK-PS-01-057	67	1.12	9.88	8.56	84.9	
K4587-02	FK-PS-01-057	68	1.12	9.65	8.41	85.5	
K4590-01	DRUM-COMP	72	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4591-01	SB-3-25-27	69	1.14	9.91	6.4	60.0	
K4591-02	SB-3-25-27	73	1.13	9.5	6.22	60.8	
K4591-03	SB-3-30-32	74	1.12	9.73	6.56	63.2	
K4591-04	SB-3-30-32	75	1.14	9.51	6.22	60.7	
K4591-06	SB-3-35-37	76	1.14	9.8	6.32	59.8	
K4591-07	SB-3-35-37	77	1.14	9.7	6.31	60.4	
K4593-01	S-3-(5-7)	78	1.14	9.85	8.62	85.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$