



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/3/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 08/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104911

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4440-06	S302-J-3.0-3.5-082019	1	1.14	9.86	8.08	79.6	
K4440-07	S302-J-4.0-4.5-082019	2	1.15	9.77	7.81	77.3	
K4440-08	S302-J-7.0-7.5-082019	3	1.12	9.58	7.61	76.7	
K4440-09	S302-J-8.0-8.5-082019	4	1.14	9.54	7.87	80.1	
K4440-10	S302-J-9.0-9.5-082019	5	1.19	9.64	8.36	84.9	
K4440-12	S303-J-2.0-2.5-082019	6	1.18	9.56	7.83	79.4	
K4440-13	S303-J-3.0-3.5-082019	7	1.15	9.96	8.13	79.2	
K4440-14	S303-J-3.0-3.5-082019-FD	8	1.16	9.97	8.16	79.5	
K4440-15	S303-J-4.0-4.5-082019	9	1.14	9.97	8.33	81.4	
K4440-16	S303-J-6.0-6.5-082019	10	1.15	9.98	8.84	87.1	
K4440-17	S303-J-7.0-7.5-082019	11	1.15	9.78	7.99	79.3	
K4440-18	S303-J-8.0-8.5-082019	12	1.15	9.94	8.45	83.0	
K4441-01	S303-J-9.0-9.5-082019	13	1.17	9.88	6.74	63.9	
K4441-03	S59-J-4.0-4.5-082019	14	1.15	9.97	8.42	82.4	
K4441-04	S59-J-5.0-5.5-082019	15	1.12	9.91	8.29	81.6	
K4441-05	S59-J-6.0-6.5-082019	16	1.18	9.97	8.36	81.7	
K4441-06	S59-J-7.0-7.5-082019	17	1.13	9.93	8.46	83.3	
K4441-07	S59-J-8.0-8.5-082019	18	1.14	9.77	8.32	83.2	
K4441-08	S59-J-9.0-9.5-082019	19	1.16	9.88	6.16	57.3	
K4441-14	S700-N-3.0-3.5-082019	20	1.1	9.84	8.27	82.0	
K4441-15	S700-N-4.0-4.5-082019	21	1.12	9.52	7.98	81.7	
K4441-16	S700-N-5.0-5.5-082019	22	1.18	9.98	8.84	87.0	
K4441-17	S700-N-6.0-6.5-082019	23	1.13	9.93	8.5	83.7	
K4441-18	S700-N-7.0-7.5-082019	24	1.14	9.97	8.48	83.1	
K4441-19	S700-N-8.0-8.5-082019	25	1.13	9.98	8.47	82.9	
K4441-20	S700-N-9.0-9.5-082019	26	1.15	9.71	8.48	85.6	
K4443-03	S701-N-2.0-2.5-082019	27	1.14	9.94	8.43	82.8	
K4443-04	S701-N-3.0-3.5-082019	28	1.14	9.61	8.18	83.1	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/3/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 08/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104911

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4443-05	S701-N-3.0-3.5-082019-FD	29	1.16	9.53	8.35	85.9	
K4443-06	S701-N-4.0-4.5-082019	30	1.1	9.61	7.95	80.5	
K4443-07	S701-N-5.0-5.5-082019	31	1.16	9.9	8.39	82.7	
K4443-08	S701-N-6.0-6.5-082019	32	1.16	9.97	8.4	82.2	
K4443-09	S701-N-7.0-7.5-082019	33	1.14	9.93	8.53	84.1	
K4443-10	S701-N-8.0-8.5-082019	34	1.16	9.86	8.5	84.4	
K4443-11	S701-N-9.0-9.5-082019	35	1.13	9.88	8.49	84.1	
K4643-01	T11-B1-B4- (0-4)	36	1.15	9.85	8.89	89.0	
K4643-02	T11-B1- (0)	37	1.11	9.61	8.7	89.3	
K4643-03	T11-B2- (0)	38	1.16	9.56	9.00	93.3	
K4643-04	T11-B3- (0)	39	1.14	9.51	8.94	93.2	
K4643-05	T11-B4- (0)	40	1.18	9.66	8.62	87.7	
K4643-06	T11-B1- (4)	41	1.13	9.53	9.07	94.5	
K4643-07	T11-B2- (4)	42	1.15	9.59	8.97	92.7	
K4643-08	T11-B3- (4)	43	1.15	9.95	9.41	93.9	
K4643-09	T11-B4- (4)	44	1.17	9.97	9.23	91.6	
K4643-10	T11-C1-C4- (0-4)	45	1.15	9.74	9.11	92.7	
K4643-11	T11-C1- (0)	46	1.18	9.85	9.16	92.0	
K4643-12	T11-C2- (0)	47	1.11	9.74	9.03	91.8	
K4643-13	T11-C3- (0)	48	1.13	9.75	8.87	89.8	
K4643-14	T11-C4- (0)	49	1.17	9.53	9.04	94.1	
K4643-15	T11-C1- (4)	50	1.18	9.79	9.22	93.4	
K4643-16	T11-C2- (4)	51	1.12	9.5	8.57	88.9	
K4643-17	T11-C3- (4)	52	1.18	9.88	9.21	92.3	
K4643-18	T11-C4- (4)	53	1.16	9.97	8.97	88.6	
K4643-19	T11-D1-D4- (0-4)	54	1.11	9.57	8.78	90.7	
K4643-20	T11-D1- (0)	55	1.13	9.81	9.01	90.8	
K4644-01	T11-D2- (0)	56	1.18	9.84	8.98	90.1	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/3/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 08/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104911

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4644-02	T11-D3- (0)	57	1.14	9.54	9.07	94.4	
K4644-03	T11-D4- (0)	58	1.14	9.91	9.13	91.1	
K4644-04	T11-D1- (4)	59	1.12	9.66	9.36	96.5	
K4644-05	T11-D2- (4)	60	1.18	9.77	9.01	91.2	
K4644-06	T11-D3- (4)	61	1.14	9.51	8.95	93.3	
K4644-07	T11-D4- (4)	62	1.18	9.71	8.99	91.6	
K4644-08	T11-E1-E4- (0-4)	63	1.14	9.73	8.43	84.9	
K4644-09	T11-E1- (0)	64	1.17	9.5	8.65	89.8	
K4644-10	T11-E2- (0)	65	1.18	9.75	8.7	87.7	
K4644-11	T11-E3- (0)	66	1.19	9.83	8.81	88.2	
K4644-12	T11-E4- (0)	67	1.15	9.96	8.73	86.0	
K4644-13	T11-E1- (4)	68	1.19	9.88	8.85	88.1	
K4644-14	T11-E2- (4)	69	1.11	9.66	8.6	87.6	
K4644-15	T11-E3- (4)	70	1.13	8.51	7.43	85.4	
K4644-16	T11-E4- (4)	71	1.1	9.81	8.45	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$