



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/4/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:43
In Date: 09/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104933

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4094-01	OILY-DEBRIS	47	1.00	2.00	2.00	100.0	OILY DEBRIS 100% SOLIDS
K4469-01	S704-N-7.0-7.5-082119	1	1.14	9.97	8.5	83.4	
K4469-02	S704-N-8.0-8.5-082119	2	1.13	9.94	8.43	82.9	
K4469-03	S704-N-9.0-9.5-082119	3	1.13	9.67	7.77	77.8	
K4470-14	S703-N-3.0-3.5-082119	4	1.14	9.77	8.47	84.9	
K4470-15	S703-N-4.0-4.5-082119	5	1.13	9.68	8.36	84.6	
K4470-16	S703-N-5.0-5.5-082119	6	1.12	9.91	7.78	75.8	
K4470-17	S703-N-6.0-6.5-082119	7	1.13	9.7	7.68	76.4	
K4470-18	S703-N-7.0-7.5-082119	8	1.13	9.73	8.22	82.4	
K4470-19	S703-N-8.0-8.5-082119	9	1.14	9.88	8.5	84.2	
K4470-20	S703-N-9.0-9.5-082119	10	1.13	9.64	8.19	83.0	
K4471-17	S704-N-3.0-3.5-082119	11	1.12	9.6	7.89	79.8	
K4471-18	S704-N-4.0-4.5-082119	12	1.14	9.9	8.35	82.3	
K4471-19	S704-N-5.0-5.5-082119	13	1.13	9.83	8.15	80.7	
K4471-20	S704-N-6.0-6.5-082119	14	1.12	9.55	8.11	82.9	
K4635-01	T-1	15	1.12	9.58	7.81	79.1	
K4635-02	T-2	16	1.13	9.77	8.61	86.6	
K4635-03	T-3	17	1.12	9.94	8.62	85.0	
K4635-04	T-4	18	1.11	9.8	8.53	85.4	
K4635-05	T-5	19	1.13	9.63	8.47	86.4	
K4639-03	DW-01	20	1.12	9.65	7.3	72.5	
K4639-04	DW-02	21	1.12	9.56	6.58	64.7	
K4639-05	DW-03	22	1.14	9.65	9.31	96.0	
K4639-06	SS-01	23	1.14	9.57	9.3	96.8	
K4653-01	287	49	1.14	9.74	8.44	84.9	
K4653-04	287	50	1.14	9.74	8.44	84.9	
K4654-01	TR-OILY-DEBRIS	48	1.00	2.00	2.00	100.0	OILY DEBRIS 100% SOLIDS
K4655-01	TP-15-S1	24	1.13	9.64	9.02	92.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/4/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:43
In Date: 09/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104933

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4655-02	TP-15-S2	25	1.15	9.95	9.11	90.5	
K4655-03	TP-15-S	26	1.13	9.54	8.79	91.1	
K4655-05	TP-15-D1	27	1.12	9.88	8.18	80.6	
K4655-06	TP-15-D2	28	1.13	9.54	8.27	84.9	
K4655-07	TP-15-D	29	1.13	9.98	8.8	86.7	
K4668-01	SU-04-090319	51	1.14	9.76	8.67	87.4	
K4668-02	SU-04-090319	52	1.13	9.88	8.05	79.1	
K4670-01	NWC-1	30	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4670-02	WC-1	31	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4670-03	DC-1	32	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4670-04	SDC-1	33	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4671-01	WC-1	34	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4672-01	DC-1	35	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4672-02	LUC-1	36	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4672-03	PJC-1	37	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4672-04	EXJC-1	38	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4673-01	DC-1	39	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4673-02	WC-1	40	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4674-01	DC-1	41	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4674-02	WC-1	42	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4674-03	EXJC-1	43	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4675-01	EDWC-1	44	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4676-01	CFM-1	45	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4676-02	DC-1	46	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$