



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 19:08

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:35

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB104986

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4470-04	S702-N-3.0-3.5-082119	1	1.15	9.56	8.44	86.7	
K4470-05	S702-N-4.0-4.5-082119	2	1.12	9.73	8.36	84.1	
K4470-06	S702-N-5.0-5.5-082119	3	1.13	9.72	8.00	80.0	
K4470-07	S702-N-6.0-6.5-082119	4	1.12	9.62	8.12	82.4	
K4470-08	S702-N-7.0-7.5-082119	5	1.12	9.64	8.35	84.9	
K4470-09	S702-N-8.0-8.5-082119	6	1.12	9.53	7.99	81.7	
K4470-10	S702-N-9.0-9.5-082119	7	1.13	9.88	8.26	81.5	
K4549-03	S700A-N-2.0-2.5-082619	8	1.13	9.77	8.69	87.5	
K4549-04	S700A-N-3.0-3.5-082619	9	1.12	9.69	8.18	82.4	
K4549-05	S700A-N-4.0-4.5-082619	10	1.12	9.63	7.87	79.3	
K4549-06	S700A-N-5.0-5.5-082619	11	1.13	9.78	7.74	76.4	
K4549-07	S700A-N-6.0-6.5-082619	12	1.12	9.85	7.87	77.3	
K4549-08	S700A-N-7.0-7.5-082619	13	1.15	9.89	7.67	74.6	
K4549-09	S700A-N-8.0-8.5-082619	14	1.12	9.66	7.78	78.0	
K4549-10	S700A-N-9.0-9.5-082619	15	1.14	9.54	7.85	79.9	
K4551-10	S400-K1-3.0-3.5-082619	16	1.13	9.8	8.52	85.2	
K4573-04	S700B-N-3.0-3.5-082719	36	1.12	9.93	8.11	79.3	
K4573-05	S700B-N-4.0-4.5-082719	37	1.14	9.53	7.7	78.2	
K4573-06	S700B-N-6.0-6.5-082719	38	1.14	9.82	8.08	80.0	
K4573-07	S700B-N-7.0-7.5-082719	39	1.13	9.67	8.26	83.5	
K4573-08	S700B-N-8.0-8.5-082719	40	1.12	9.59	7.97	80.9	
K4573-09	S700B-N-9.0-9.5-082719	41	1.12	9.83	8.94	89.8	
K4573-15	S702A-N-3.0-3.5-082719	42	1.12	9.54	7.1	71.0	
K4573-16	S702A-N-4.0-4.5-082719	43	1.14	9.58	7.9	80.1	
K4573-17	S702A-N-6.0-6.5-082719	44	1.14	9.93	8.39	82.5	
K4573-18	S702A-N-7.0-7.5-082719	45	1.15	9.74	8.06	80.4	
K4574-10	S702A-N-8.0-8.5-082719	46	1.14	9.71	8.11	81.3	
K4574-11	S702A-N-9.0-9.5-082719	47	1.13	9.82	8.38	83.4	



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 19:08

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:35

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB104986

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4575-16	S703A-N-3.0-3.5-082719	48	1.14	9.92	8.54	84.3	
K4575-17	S703A-N-4.0-4.5-082719	49	1.13	9.63	7.93	80.0	
K4575-18	S703A-N-6.0-6.5-082719	50	1.14	9.9	8.11	79.6	
K4575-19	S703A-N-7.0-7.5-082719	51	1.14	9.84	8.4	83.4	
K4575-20	S703A-N-8.0-8.5-082719	52	1.15	9.6	8.14	82.7	
K4575-21	S703A-N-9.0-9.5-082719	53	1.12	9.7	8.08	81.1	
K4658-01	HR-06-090619	54	1.13	9.62	8.91	91.6	
K4658-02	HR-06-090619	55	1.13	9.6	8.65	88.8	
K4658-03	HR-06-090619-1	56	1.14	9.65	8.53	86.8	
K4658-04	HR-06-090619-2	57	1.12	9.91	8.77	87.0	
K4658-05	HR-06-090619-3	58	1.12	9.7	8.62	87.4	
K4662-01	CL-01-090619	59	1.13	9.66	9.11	93.6	
K4662-02	CL-01-090619	60	1.12	9.87	9.3	93.5	
K4662-03	CL-01-090619-1	61	1.12	9.65	8.83	90.4	
K4662-04	CL-01-090619-2	62	1.13	9.64	9.09	93.5	
K4662-05	CL-01-090619-3	63	1.12	9.83	9.28	93.7	
K4662-06	CL-01-090619-4	64	1.14	9.73	8.88	90.1	
K4662-07	CL-01-090619-5	65	1.12	9.66	8.94	91.6	
K4664-01	OR-03-090619	66	1.11	9.78	9.14	92.6	
K4664-02	OR-03-090619	67	1.12	9.76	9.38	95.6	
K4664-03	OR-03-090619-1	68	1.12	9.75	9.15	93.0	
K4664-04	OR-03-090619-2	69	1.12	9.88	9.11	91.2	
K4719-01	T4-1-4- (0-1.5)	17	1.12	9.75	8.39	84.2	
K4719-02	T4-1- (0)	18	1.12	9.78	8.86	89.4	
K4719-03	T4-1- (1.5)	19	1.13	9.97	9.32	92.6	
K4719-04	T4-2- (0)	20	1.13	9.97	8.68	85.4	
K4719-05	T4-2- (1.5)	21	1.12	9.86	8.56	85.1	
K4719-06	T4-3- (0)	22	1.11	9.29	8.1	85.5	



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 19:08

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:35

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB104986

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4719-07	T4-3-(1.5)	23	1.12	8.25	6.93	81.5	
K4719-08	T4-4-(0)	24	1.13	9.88	8.46	83.8	
K4719-09	T4-4-(1.5)	25	1.12	9.74	8.06	80.5	
K4719-10	T4-1-4-(0-1.5)	26	1.12	9.6	9.04	93.4	
K4719-11	T4-2-(0)	27	1.12	9.86	8.96	89.7	
K4719-12	T4-2-(1.5)	28	1.14	9.54	8.65	89.4	
K4719-13	T4-3-(0)	29	1.14	9.8	9.41	95.5	
K4719-14	T4-3-(1.5)	30	1.12	9.94	9.45	94.4	
K4719-15	T4-4-(0)	31	1.14	9.8	9.05	91.3	
K4719-16	T4-4-(1.5)	32	1.14	9.6	8.82	90.8	
K4721-01	S700-A2-N-0.0-0.5-090519	33	1.12	9.86	6.89	66.0	
K4721-11	S702-B1-N-0.0-0.5-090519	34	1.12	9.6	9.1	94.1	
K4722-01	S703-A1-N-0.0-0.5-090519	35	1.12	9.55	8.58	88.5	
K4753-01	DRUM-1-OIL	70	1.00	2.00	2.00	100.0	OIL SAMPLE
K4753-03	DRUM-2-RAGS	71	1.00	2.00	2.00	100.0	100 % SOLIDS
K4755-01	RO-1-3	72	1.14	9.88	8.09	79.5	
K4755-02	RO-1-3	73	1.14	9.65	8.16	82.5	
K4757-01	A-0-4	74	1.14	9.77	8.42	84.4	
K4757-02	A-0-4	75	1.13	6.61	5.59	81.4	
K4757-03	A-4-8	76	1.14	9.7	7.64	75.9	
K4757-04	A-4-8	77	1.15	9.75	8.13	81.2	
K4757-05	A1	78	1.15	9.67	7.77	77.7	
K4757-06	A2	79	1.15	9.68	7.94	79.6	
K4757-07	A3	80	1.13	9.91	7.96	77.8	
K4757-08	A4	81	1.12	9.58	7.81	79.1	
K4757-09	A5	82	1.13	9.67	8.07	81.3	
K4757-10	A6	83	1.12	9.74	8.09	80.9	



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 19:08

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:35

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB104986

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4757-11	A7	84	1.14	9.82	8.03	79.4	
K4757-12	A8	85	1.13	9.55	7.89	80.3	
K4757-13	A9	86	1.12	9.89	8.29	81.8	
K4757-14	A10	87	1.12	9.91	8.16	80.1	
K4757-15	A11	88	1.13	9.93	8.54	84.2	
K4757-16	A12	89	1.14	9.87	8.06	79.3	
K4757-17	A13	90	1.13	9.6	8.11	82.4	
K4757-18	A14	91	1.12	9.72	8.16	81.9	
K4757-19	A15	92	1.12	9.7	8.17	82.2	
K4757-20	A16	93	1.12	9.66	8.00	80.6	
K4758-01	B-0-4	94	1.13	9.73	8.5	85.7	
K4758-02	B-0-4	95	1.12	5.13	4.6	86.8	
K4758-03	B-4-8	96	1.14	9.94	9.28	92.5	
K4758-04	B-4-8	97	1.15	9.82	9.19	92.7	
K4758-05	B1	98	1.14	9.6	8.44	86.3	
K4758-06	B2	99	1.14	9.58	8.98	92.9	
K4758-07	B3	100	1.14	9.64	8.86	90.8	
K4758-08	B4	101	1.15	9.77	8.6	86.4	
K4758-09	B5	102	1.15	9.65	8.87	90.8	
K4758-10	B6	103	1.12	9.71	8.49	85.8	
K4758-11	B7	104	1.14	9.84	8.94	89.7	
K4758-12	B8	105	1.13	9.82	8.96	90.1	
K4758-13	B9	106	1.14	9.85	8.97	89.9	
K4758-14	B10	107	1.12	9.79	8.8	88.6	
K4758-15	B11	108	1.16	9.58	8.36	85.5	
K4758-16	B12	109	1.15	9.88	9.11	91.2	
K4758-17	B13	110	1.13	9.98	9.12	90.3	
K4758-18	B14	111	1.14	9.72	9.00	91.6	



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 19:08

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:35

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB104986

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4758-19	B15	112	1.14	9.81	8.9	89.5	
K4758-20	B16	113	1.12	9.76	9.14	92.8	
K4759-01	C-0-4	114	1.12	9.86	9.05	90.7	
K4759-02	C-0-4	115	1.14	9.95	9.22	91.7	
K4759-03	C-4-8	116	1.13	9.76	8.8	88.9	
K4759-04	C-4-8	117	1.12	9.98	8.83	87.0	
K4759-05	C1	118	1.13	9.74	8.56	86.3	
K4759-06	C2	119	1.13	9.94	8.97	89.0	
K4759-07	C3	120	1.13	9.58	8.72	89.8	
K4759-08	C4	121	1.14	9.56	8.75	90.4	
K4759-09	C5	122	1.12	9.67	8.81	89.9	
K4759-10	C6	123	1.11	9.76	8.93	90.4	
K4759-11	C7	124	1.13	9.92	8.7	86.1	
K4759-12	C8	125	1.13	9.95	9.1	90.4	
K4759-13	C9	126	1.13	9.55	8.76	90.6	
K4759-14	C10	127	1.12	9.81	8.96	90.2	
K4759-15	C11	128	1.13	9.93	9.06	90.1	
K4759-16	C12	129	1.12	9.83	9.05	91.0	
K4759-17	C13	130	1.12	9.9	9.02	90.0	
K4759-18	C14	131	1.13	9.95	9.06	89.9	
K4759-19	C15	132	1.12	9.93	8.97	89.1	
K4759-20	C16	133	1.12	9.82	8.75	87.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$