



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 107

Time IN: 15:00

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 07:20

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB105000

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4741-01	MH0ZR5	1	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-02	MH1148	2	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-03	MH1148D	3	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-04	MH1148S	4	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-05	MH11T9	5	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-06	MH11W0	6	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-07	MH11W4	7	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-08	MH11W5	8	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-09	MH11W6	9	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-10	MH11W7	10	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-11	MH11W8	11	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-12	MH11X0	12	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-13	MH11X2	13	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-14	MH11X3	14	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-15	MH11X4	15	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-16	MH11X5	16	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-17	MH11X6	17	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-18	MH11X7	18	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-19	MH11Y1	19	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-20	MH11Y2	20	1.00	2.00	2.00	100.0	100% SOLIDS
K4741-21	MH11Y3	21	1.00	2.00	2.00	100.0	100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$