



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/10/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:05
In Date: 09/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:26
Out Date: 09/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105001

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4743-01	MH11Y4	1	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-02	MH11Y5	2	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-03	MH11Y6	3	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-04	MH11Y7	4	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-05	MH11Y8	5	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-06	MH11Y8D	6	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-07	MH11Y8S	7	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-08	MH12D5	8	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-09	MH12D6	9	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-10	MH12D7	10	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-11	MH12D8	11	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-12	MH12E0	12	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-13	MH12E1	13	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-14	MH12E2	14	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-15	MH12E3	15	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-16	MH12E5	16	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-17	MH12E6	17	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-18	MH12E7	18	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-19	MH12E8	19	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-20	MH12E9	20	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-21	MH12F0	21	1.00	2.00	2.00	100.0	100% SOLIDS
K4743-22	MH12F1	22	1.00	2.00	2.00	100.0	100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$