



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 9/9/2019

OVENTEMP IN Celsius(°C): 107

Time IN: 15:20

In Date: 09/06/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 07:30

Out Date: 09/07/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB105003

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4748-01	MH12M6	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-02	MH12M7	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-03	MH12M8	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-04	MH12M9	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-05	MH12N0	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-06	MH12N1	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-07	MH12N2	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-08	MH12N3	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-09	MH12N4	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-10	MH12N4D	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-11	MH12N4S	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-12	MH12N5	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-13	MH12N6	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-14	MH12N7	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-15	MH12N8	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-16	MH12N9	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-17	MH12P0	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-18	MH12P1	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-19	MH12P2	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-20	MH12P3	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4748-21	MH12P4	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$