



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/10/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:32
In Date: 09/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:45
Out Date: 09/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105041

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4697-13	S400B-K1-2.0-2.5-090419	1	1.15	9.6	8.63	88.5	
K4698-14	S400D-K1-2.0-2.5-090419	2	1.14	9.8	8.35	83.3	
K4738-01	FK-BPM-01-H02-H03-H02A-H03A	3	1.13	9.7	9.02	92.1	
K4738-02	FK-BPM-01-H02-H03-H02A-H03A	4	1.14	9.73	9.04	92.0	
K4762-10	STOCK-PILE	5	1.14	9.83	8.65	86.4	
K4764-01	ROLL-OFF-VOC-1	6	1.11	9.7	8.13	81.7	
K4764-02	ROLL-OFF-VOC-2	7	1.12	9.82	8.93	89.8	
K4764-03	ROLL-OFF-VOC-3	8	1.14	9.83	9.37	94.7	
K4764-04	ROLL-OFF-COMPOSITE-1	9	1.13	9.82	7.94	78.4	
K4764-05	ROLL-OFF-COMPOSITE-2	10	1.14	9.88	8.93	89.1	
K4764-06	ROLL-OFF-COMPOSITE-3	11	1.15	9.74	8.89	90.1	
K4764-07	ROLL-OFF-TPH-1	12	1.14	9.67	8.83	90.2	
K4764-08	ROLL-OFF-TPH-2	13	1.14	9.66	8.9	91.1	
K4764-09	ROLL-OFF-TPH-3	14	1.12	9.7	8.8	89.5	
K4764-14	ROLL-OFF-COMPOSITE-1	15	1.11	9.7	8.13	81.7	
K4764-15	ROLL-OFF-COMPOSITE-2	16	1.12	9.82	8.93	89.8	
K4764-16	ROLL-OFF-COMPOSITE-3	17	1.14	9.83	9.37	94.7	
K4766-01	390903706	18	1.14	9.88	1.53	4.5	
K4767-01	390903773	19	1.15	9.51	3.8	31.7	
K4768-01	390905716	20	1.15	9.96	2.59	16.3	
K4786-01	TP1	21	1.14	9.58	9.02	93.4	
K4786-02	TP1-B	22	1.15	9.7	9.09	92.9	
K4786-03	TP2	23	1.15	9.86	9.24	92.9	
K4786-04	TP2-B	24	1.14	9.82	9.15	92.3	
K4786-05	TP3	25	1.15	9.8	9.27	93.9	
K4786-07	TP4	26	1.14	9.64	9.26	95.5	
K4786-08	TP-4	27	1.16	9.72	9.28	94.9	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/10/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:32
In Date: 09/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:45
Out Date: 09/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105041

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4787-01	TP-3-S(1-4)	28	1.15	9.86	9.19	92.3	
K4787-02	TP-3-S1(1-1.5)	29	1.15	9.89	9.28	93.0	
K4787-03	TP-3-S2(3.5-4)	30	1.14	9.73	8.94	90.8	
K4787-05	TP-5-S(1-4)	31	1.16	9.77	8.88	89.7	
K4787-06	TP-5-S1(1.5-2)	32	1.16	9.91	8.95	89.0	
K4787-07	TP-5-S2(3.5-4)	33	1.15	9.94	9.1	90.4	
K4787-09	TP-5-D(4-10)	34	1.14	9.87	8.82	88.0	
K4787-10	TP-5-D1-(5-6)	35	1.14	9.56	8.66	89.3	
K4787-11	TP-5-D2(9.5-10)	36	1.14	9.98	8.62	84.6	
K4789-01	RT3465	37	1.12	9.9	9.45	94.9	
K4789-02	RT3465	38	1.12	9.77	9.51	97.0	
K4789-04	RT5124	39	1.13	9.6	9.07	93.7	
K4789-05	RT5124	40	1.14	9.92	8.79	87.1	
K4789-07	RT3407	41	1.14	9.7	9.15	93.6	
K4789-08	RT3407	42	1.13	9.76	7.79	77.2	
K4789-10	RT3860	43	1.14	9.9	8.22	80.8	
K4789-12	RT3860	44	1.15	9.92	9.63	96.7	
K4789-14	RBR200029	45	1.14	9.87	9.22	92.6	
K4789-15	RBR200029	46	1.12	9.85	9.42	95.1	
K4790-04	DRUM-2	47	1.14	9.55	9.39	98.1	
K4790-07	DRUM-3	48	1.14	9.55	9.11	94.8	
K4791-01	RBR200030	49	1.14	9.55	8.55	88.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$