



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/11/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:30
In Date: 09/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105059

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K3631-01	WALL-PANEL	1	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K3631-03	WALL-BASE	2	1.00	2.00	2.00	100.0	STONE SAMPLE 100% SOLIDS
K4551-11	S400-K1-4.0-4.5-082619	3	1.16	9.92	8.62	85.2	
K4721-02	S700-A2-N-1.0-1.5-090519	4	1.14	9.95	8.83	87.3	
K4804-01	VB-16132	5	1.14	9.98	6.49	60.5	SLUDGE SAMPLE
K4805-01	RED-FILL-DIRT	6	1.14	9.85	9.48	95.8	
K4806-01	DRUM-COMP	7	1.13	9.59	5.45	51.1	
K4814-01	WC-1	8	1.15	9.92	8.72	86.3	
K4814-02	WC-2	9	1.13	9.7	8.44	85.3	
K4814-03	WC-3	10	1.13	9.7	8.58	86.9	
K4819-01	TP-2-S1 (1-1.5)	11	1.15	9.86	9.3	93.6	
K4819-02	TP-2-S2 (3.5-4)	12	1.15	9.78	7.54	74.0	
K4819-03	TP-2-S (1-4)	13	1.13	9.98	9.02	89.2	
K4819-05	TP-7-S1 (2-2.5)	14	1.12	9.66	8.71	88.9	
K4819-06	TP-7-S2 (3.5-4)	15	1.12	9.85	8.9	89.1	
K4819-07	TP-7-S (2-4)	16	1.14	9.96	9.22	91.6	
K4819-09	TP-7-D1 (5-5.5)	17	1.14	9.78	9.35	95.0	
K4819-10	TP-7-D2 (9.5-10)	18	1.14	9.98	9.34	92.8	
K4819-11	TP-7-D (4-10)	19	1.13	9.58	8.79	90.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$