



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/12/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 09/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 09/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105071

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4779-01	C0AF7	1	1.13	9.57	8.81	91.0	
K4779-02	C0AF8	2	1.13	9.9	8.97	89.4	
K4779-03	C0AG2	3	1.15	9.7	8.53	86.3	
K4779-04	C0AG3	4	1.13	9.67	8.7	88.6	
K4779-05	C0AG4	5	1.14	9.96	7.86	76.2	
K4779-06	C0AG5	6	1.18	9.69	7.68	76.4	
K4779-07	C0AG6	7	1.15	9.95	8.33	81.6	
K4779-08	C0AG7	8	1.18	9.74	8.23	82.4	
K4779-09	C0AG8	9	1.17	9.59	7.79	78.6	
K4779-10	C0AG9	10	1.15	9.88	8.32	82.1	
K4779-11	C0AH9	11	1.16	9.76	8.77	88.5	
K4779-12	C0AJ0	12	1.15	9.96	8.78	86.6	
K4779-13	C0AJ1	13	1.16	9.65	8.69	88.7	
K4779-14	C0AJ4	14	1.14	9.71	8.22	82.6	
K4779-15	C0AJ5	15	1.14	9.64	8.39	85.3	
K4779-16	C0AJ6	16	1.13	9.97	8.86	87.4	
K4779-17	C0AJ7	17	1.1	9.79	7.96	78.9	
K4779-18	C0AJ8	18	1.18	9.67	8.05	80.9	
K4779-19	C0AJ9	19	1.14	9.98	8.58	84.2	
K4779-20	C0AK0	20	1.14	9.9	7.21	69.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$