



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:20
In Date: 09/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105096

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4663-04	OK-01-091119	37	1.15	9.96	9.7	97.0	
K4663-05	OK-01-091119	38	1.14	9.87	9.62	97.1	
K4666-01	PL-03-091219-A	67	1.13	9.84	9.39	94.8	
K4666-03	PL-03-091219-B	68	1.13	9.75	9.56	97.8	
K4666-05	PL-04-091219-1	69	1.15	9.91	9.22	92.1	
K4666-06	PL-04-091219-2	70	1.12	9.55	9.15	95.3	
K4666-07	PL-04-091219-3	71	1.15	9.7	9.04	92.3	
K4666-08	PL-04-091219-4	72	1.16	9.57	9.51	99.3	
K4666-09	PL-04-091219-5	73	1.12	9.95	9.13	90.7	
K4666-10	PL-04-091219-6	74	1.15	9.55	8.99	93.3	
K4666-11	PL-04-091219-7	75	1.14	9.84	8.53	84.9	
K4668-04	SU-02-091219	39	1.12	9.72	9.06	92.3	
K4668-05	SU-02-091219-1	40	1.14	9.88	9.67	97.6	
K4668-06	SU-02-091219-2	41	1.14	9.88	9.67	97.6	
K4846-01	SB-19-15.5-16.0	1	1.15	9.74	9.44	96.5	
K4846-02	SB-19-COMP	2	1.14	9.87	9.34	93.9	
K4846-03	SB-21-13.5-14.0	3	1.13	9.61	8.2	83.4	
K4846-04	SB-21-COMP	4	1.13	9.69	8.31	83.9	
K4846-05	SB-22-7.5-8.0	5	1.13	9.86	9.09	91.2	
K4846-06	SB-22-COMP	6	1.14	9.86	9.08	91.1	
K4846-07	SB-23-5.5-6.0	7	1.12	9.67	8.77	89.5	
K4846-08	SB-23-COMP	8	1.12	9.7	8.82	89.7	
K4846-09	SB-25-7.5-8.0	9	1.12	9.8	8.48	84.8	
K4846-10	SB-25-COMP	10	1.12	9.57	8.5	87.3	
K4847-01	SB-27-7.5-8.0	11	1.12	9.71	8.94	91.0	
K4847-02	SB-27-COMP	12	1.13	9.6	8.96	92.4	
K4847-03	SB-28-5.5-6.0	13	1.13	9.75	8.59	86.5	
K4847-04	SB-28-COMP	14	1.13	9.78	8.87	89.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:20
In Date: 09/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105096

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4847-05	SB-29-12.5-13.0	15	1.12	9.92	8.96	89.1	
K4847-06	SB-29-COMP	16	1.12	9.81	8.98	90.4	
K4847-07	SB-34-13.5-14.0	17	1.14	9.7	9.06	92.5	
K4847-08	SB-34-COMP	18	1.14	9.65	8.92	91.4	
K4847-09	SB-36-10.5-11.0	19	1.14	9.5	8.86	92.3	
K4847-10	SB-36-COMP	20	1.13	9.64	8.9	91.3	
K4847-11	SB-35-5.5-6.0	21	1.13	9.53	9.03	94.0	
K4847-12	SB-35-COMP	22	1.14	9.86	9.37	94.4	
K4847-13	SB-37-10.5-11.0	23	1.12	9.7	8.71	88.5	
K4847-14	SB-37-COMP	24	1.12	9.8	8.93	90.0	
K4851-01	SB-3C- (9)	25	1.13	8.43	6.88	78.8	
K4851-02	K4851-01MS	26	1.13	8.43	6.88	78.8	
K4851-03	K4851-01MSD	27	1.13	8.43	6.88	78.8	
K4851-04	SB-3C- (9-10)	28	1.14	9.68	7.62	75.9	
K4851-06	SB-4C- (9)	29	1.14	9.71	8.49	85.8	
K4851-07	SB-4C- (9-10)	30	1.13	9.57	8.8	90.9	
K4851-09	SB-3D- (9)	31	1.13	9.76	8.12	81.0	
K4851-10	SB-3D- (9-10)	32	1.13	9.82	8.06	79.7	
K4861-01	TP-11-S1 (1.5-2)	76	1.14	9.56	8.91	92.3	
K4861-02	TP-11-S2 (3.5-4)	77	1.14	9.64	9.11	93.8	
K4861-03	TP-11-S (1.5-4)	78	1.14	9.96	8.57	84.2	
K4861-05	TP-11-D1 (4-4.5)	79	1.15	9.73	8.94	90.8	
K4861-06	TP-11-D2 (6-6.5)	80	1.15	9.6	9.22	95.5	
K4861-07	TP-11-D (4-6.5)	81	1.14	9.85	8.06	79.4	
K4861-09	TP-1-S1 (1.5-2)	82	1.14	9.52	8.44	87.1	
K4861-10	TP-1-S2 (3.5-4)	83	1.14	9.65	8.36	84.8	
K4861-11	TP-1-S (1.5-4)	84	1.14	9.73	8.69	87.9	
K4862-02	CSA-4-1	42	1.13	9.92	9.67	97.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:20
In Date: 09/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105096

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4862-03	CSA-4-2	43	1.14	9.78	9.55	97.3	
K4862-04	CSA4-3	44	1.13	9.51	9.24	96.8	
K4862-05	CSA4-4	45	1.11	9.7	9.41	96.6	
K4862-06	CSA4-4	46	1.15	9.73	9.58	98.3	
K4862-07	CSA4-4-.05	47	1.13	9.65	9.53	98.6	
K4862-08	CSA4-4-1	48	1.15	9.54	9.42	98.6	
K4862-09	CSA4-1.5	49	1.15	9.51	9.41	98.8	
K4862-10	CSA4-5	50	1.13	9.51	8.91	92.8	
K4863-01	CSA-4-5	51	1.13	9.65	9.01	92.5	
K4863-02	CSA-4-5-1	52	1.16	9.68	9.22	94.6	
K4863-03	CSA-4-5-2	53	1.14	9.52	8.72	90.5	
K4863-04	CSA-4-6	54	1.13	9.8	9.31	94.3	
K4863-05	CSA-4-6	55	1.14	9.6	9.1	94.1	
K4863-06	CSA-4-6-0.5	56	1.13	9.78	9.56	97.5	
K4863-07	CSA-4-6-1	57	1.14	9.54	9.11	94.9	
K4863-08	CSA-4-6-1.5	58	1.13	9.53	9.01	93.8	
K4863-09	CSA-4-7	59	1.14	9.58	8.68	89.3	
K4863-10	CSA-4-7	60	1.14	9.61	8.96	92.3	
K4863-11	CSA-4-7-1	61	1.12	9.9	9.31	93.3	
K4863-12	CSA-4-7-2	62	1.15	9.74	8.76	88.6	
K4863-13	CSA-2-5	63	1.14	9.58	8.93	92.3	
K4863-14	CSA-2-5	64	1.14	9.65	8.85	90.6	
K4863-15	CSA-2-5-2	65	1.14	9.54	8.85	91.8	
K4863-16	CSA-2-5-4	66	1.12	9.72	9.08	92.6	
K4864-01	SPUD-1-091119	33	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4864-02	SPUD-2-091119	34	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4864-03	SPUD-3-091119	35	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS
K4864-04	SPUD-4-091119	36	1.00	2.00	2.00	100.0	WIPE SAMPLE 100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:20
In Date: 09/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105096

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4866-04	ROLL-OFF-1	85	1.14	9.88	8.31	82.0	
K4866-05	ROLL-OFF-2	86	1.14	9.56	9.22	96.0	
K4873-01	VOC-28	87	1.14	5.03	4.64	90.0	
K4873-02	COMP-28	88	1.14	9.56	8.68	89.5	
K4873-03	VOC-29	89	1.15	4.65	4.47	94.9	
K4873-04	COMP-29	90	1.14	9.79	9.22	93.4	
K4873-05	VOC-30	91	1.14	4.29	4.19	96.8	
K4873-06	COMP-30	92	1.14	9.89	9.33	93.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$