



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:20
In Date: 09/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105108

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4845-01	C0AM3	1	1.12	9.66	8.06	81.3	
K4845-02	C0AA1	2	1.14	9.89	7.85	76.7	
K4845-03	C0AM5	3	1.12	9.5	7.84	80.2	
K4845-04	C0AM6	4	1.12	9.86	8.37	83.0	
K4845-05	C0AM7	5	1.14	9.81	8.39	83.6	
K4845-06	C0AM8	6	1.12	9.51	7.92	81.0	
K4845-07	C0AM9	7	1.12	9.91	8.21	80.7	
K4845-08	C0AN0	8	1.12	9.7	8.62	87.4	
K4845-09	C0AN1	9	1.12	9.6	7.59	76.3	
K4845-10	C0AN2	10	1.13	9.74	8.35	83.9	
K4845-11	C0AN3	11	1.12	9.64	7.41	73.8	
K4845-12	C0AN4	12	1.11	9.8	8.24	82.0	
K4845-13	C0AN5	13	1.12	9.6	7.78	78.5	
K4845-14	C0AN6	14	1.12	9.64	8.31	84.4	
K4845-15	C0AN7	15	1.13	9.65	8.37	85.0	
K4845-16	C0AN9	16	1.13	9.98	8.6	84.4	
K4845-17	C0AN9MS	17	1.13	9.98	8.6	84.4	
K4845-18	C0AN9MSD	18	1.13	9.98	8.6	84.4	
K4845-19	C0AP1	19	1.14	9.6	8.33	85.0	
K4845-20	C0AP2	20	1.12	9.98	7.9	76.5	
K4845-21	C0AP3	21	1.14	9.64	8.32	84.5	
K4845-22	C0AP4	22	1.15	9.7	8.36	84.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$