



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105112

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4596-02	S403-K1-2.0-2.5-082819	126	1.14	9.72	8.39	84.5	
K4661-04	HD-02-091319	97	1.14	9.64	8.82	90.4	
K4661-05	HD-02-091319	98	1.12	9.6	9.15	94.7	
K4664-06	OR-02-091319	99	1.16	9.92	9.26	92.5	
K4664-07	OR-02-091319-1	100	1.16	9.88	9.09	90.9	
K4664-08	OR-02-091319-2	101	1.14	9.96	9.03	89.5	
K4853-01	RI14-7- (24-30)	1	1.16	9.56	8.26	84.5	
K4853-02	RI14-4- (24-30)	2	1.12	9.67	7.98	80.2	
K4853-03	RI2-2- (30-36)	3	1.16	9.52	8.06	82.5	
K4853-04	AOC14-02-B4- (24-30)	4	1.15	9.55	8.43	86.7	
K4853-05	RI2-5- (30-36)	5	1.13	9.7	8.93	91.0	
K4853-06	RI14-3- (24-30)	6	1.14	9.6	8.12	82.5	
K4853-07	RI12-8- (30-36)	7	1.14	9.6	8.14	82.7	
K4853-08	RI2-1- (30-36)	8	1.14	9.73	7.86	78.2	
K4853-09	RI14-9- (24-30)	9	1.15	9.91	8.25	81.1	
K4853-10	RI14-06- (24-30)	10	1.15	9.66	8.14	82.1	
K4853-11	RI14-5- (30-36)	11	1.13	9.73	8.54	86.2	
K4853-12	RI14-5- (24-30)	12	1.15	9.72	8.39	84.5	
K4853-13	RI2-3- (30-36)	13	1.15	9.62	8.28	84.2	
K4853-14	RI14-08- (24-30)	14	1.14	9.52	8.27	85.1	
K4853-15	RI2-6- (30-36)	15	1.14	9.63	8.86	90.9	
K4853-16	RI2-7- (30-36)	16	1.13	9.52	7.51	76.0	
K4854-01	RI2-5- (24-30)	17	1.13	9.88	9.22	92.5	
K4854-02	RI14-9- (0-6)	18	1.14	9.81	9.51	96.5	
K4854-03	RI14-9- (6-12)	19	1.12	9.74	8.55	86.2	
K4854-04	RI14-9- (12-18)	20	1.14	9.54	9.00	93.6	
K4854-05	RI14-9- (18-24)	21	1.14	9.72	8.5	85.8	
K4854-06	RI14-6- (0-6)	22	1.13	9.8	9.00	90.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105112

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4854-07	RI14-06-(6-12)	23	1.12	9.65	8.55	87.1	
K4854-08	RI14-06-(12-18)	24	1.16	9.56	8.46	86.9	
K4854-09	RI14-06-(18-24)	25	1.14	9.52	8.15	83.7	
K4854-10	RI2-1-(24-30)	26	1.12	9.67	8.12	81.9	
K4854-11	RI2-8-(0-6)	27	1.15	9.63	7.77	78.1	
K4854-12	RI2-8-(6-12)	28	1.14	9.67	8.76	89.3	
K4854-13	RI2-8-(12-18)	29	1.15	9.51	8.47	87.6	
K4854-14	RI2-8-(18-24)	30	1.15	9.53	8.29	85.2	
K4854-15	RI2-8-(24-30)	31	1.13	9.78	8.44	84.5	
K4854-16	RI2-6-(0-6)	32	1.13	9.84	9.3	93.8	
K4854-17	RI2-6-(6-12)	33	1.16	9.72	9.14	93.2	
K4854-18	RI2-6-(12-18)	34	1.12	9.78	9.09	92.0	
K4854-19	RI2-6-(18-24)	35	1.15	9.63	9.00	92.6	
K4854-20	RI2-6-(24-30)	36	1.13	9.74	8.81	89.2	
K4855-01	RI2-7-(0-6)	37	1.15	9.56	8.21	83.9	
K4855-02	RI2-7-(6-12)	38	1.12	9.72	8.73	88.5	
K4855-03	RI2-7-(12-18)	39	1.12	9.7	8.42	85.1	
K4855-04	RI2-7-(18-24)	40	1.13	9.83	8.48	84.5	
K4855-05	RI2-7-(24-30)	41	1.13	9.58	8.26	84.4	
K4855-06	RI2-5-(0-6)	42	1.14	9.76	9.16	93.0	
K4855-07	RI2-5-(6-12)	43	1.16	9.66	9.06	92.9	
K4855-08	RI2-5-(12-18)	44	1.14	9.8	9.15	92.5	
K4855-09	RI2-5-(18-24)	45	1.14	9.82	9.22	93.1	
K4855-10	RI14-07-(0-6)	46	1.14	9.83	9.17	92.4	
K4855-11	RI14-07-(6-12)	47	1.12	9.68	8.94	91.4	
K4855-12	RI14-07-(12-18)	48	1.14	9.81	8.39	83.6	
K4855-13	RI14-07-(18-24)	49	1.13	9.83	8.66	86.6	
K4855-14	RI14-4-(18-24)	50	1.13	9.8	8.61	86.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105112

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4855-15	RI14-08-(18-24)	51	1.14	9.62	7.63	76.5	
K4855-16	RI14-08-(12-18)	52	1.14	9.62	8.59	87.9	
K4855-17	RI14-3-(18-24)	53	1.13	9.75	8.53	85.8	
K4855-18	AOC14-0804-(12-18)	54	1.14	9.7	8.58	86.9	
K4855-19	AOC14-0804-(18-24)	55	1.13	9.54	8.43	86.8	
K4856-01	RI14-08-(0-6)	56	1.12	9.8	9.08	91.7	
K4856-02	RI14-08-(6-12)	57	1.13	9.65	8.69	88.7	
K4856-03	RI14-05-(18-24)	58	1.14	9.74	8.32	83.5	
K4856-04	RI2-2-(24-30)	59	1.13	9.62	8.43	86.0	
K4856-05	RI2-3-(24-30)	60	1.14	9.54	8.3	85.2	
K4856-06	DUP-1	61	1.13	9.96	9.22	91.6	
K4856-07	DUP-2	62	1.13	9.77	9.44	96.2	
K4856-08	DUP-3	63	1.15	9.63	9.03	92.9	
K4860-01	IDW-SOIL-091119	64	1.14	9.53	7.78	79.1	
K4877-02	CSA2-2	65	1.13	9.8	8.94	90.1	
K4877-03	CSA2-2-1	66	1.15	9.87	8.63	85.8	
K4877-04	CSA2-2-2	67	1.14	9.52	8.7	90.2	
K4877-05	CSA2-2-3	68	1.15	9.9	8.81	87.5	
K4877-06	CSA2-3	69	1.12	9.85	8.57	85.3	
K4877-07	CSA2-3-2	70	1.15	9.64	8.92	91.5	
K4877-08	CSA2-3-4	71	1.13	9.65	8.79	89.9	
K4877-09	CSA2-4	72	1.14	9.66	9.05	92.8	
K4877-10	CSA2-4-1	73	1.14	9.68	9.00	92.0	
K4878-01	CSA-2-4-2	74	1.15	9.87	9.24	92.8	
K4878-02	CSA-2-4-3	75	1.14	9.77	9.09	92.1	
K4878-03	CSA-2-6	76	1.14	9.83	9.03	90.8	
K4878-04	CSA-2-6-2	77	1.14	9.88	9.33	93.7	
K4878-05	CSA-2-6-4	78	1.14	9.56	8.68	89.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105112

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4878-06	CSA-2-7	79	1.12	9.73	8.9	90.4	
K4878-07	CSA-2-7-2	80	1.12	9.54	8.82	91.4	
K4878-08	CSA-2-7-4	81	1.15	9.66	8.6	87.5	
K4878-09	CSA-2-8	82	1.13	9.58	8.55	87.8	
K4878-10	CSA-2-8-1	83	1.15	9.67	9.00	92.1	
K4878-11	CSA-2-8-2	84	1.14	9.84	8.61	85.9	
K4878-12	CSA-2-8-3	85	1.15	9.68	8.36	84.5	
K4878-13	CSA-2-9	86	1.13	9.6	8.49	86.9	
K4878-14	CSA-2-9-2	87	1.12	9.75	8.44	84.8	
K4878-15	CSA-2-9-4	88	1.14	9.55	8.41	86.4	
K4878-16	TPA-9	89	1.15	9.61	8.83	90.8	
K4878-17	TPA-9	90	1.15	9.88	8.97	89.6	
K4878-26	CSA-2-1	91	1.15	9.8	9.12	92.1	
K4878-27	CSA-2-1-2	92	1.14	9.54	8.78	91.0	
K4878-28	CSA-4-1-4	93	1.15	9.68	8.99	91.9	
K4884-01	6INCH-PIT	102	1.13	9.88	9.19	92.1	
K4885-01	RT-1868-RBR-200035	103	1.14	9.6	8.02	81.3	
K4886-01	TP-13-S1 (1.5-2)	104	1.14	9.71	8.82	89.6	
K4886-02	TP-13-S2 (3.5-4)	105	1.13	9.86	9.06	90.8	
K4886-03	TP-13-S (1.5-4)	106	1.13	9.74	9.04	91.9	
K4888-01	979-S-1- (18-19)	94	1.13	9.8	9.14	92.4	
K4888-02	980-B-1- (24)	95	1.13	9.57	8.21	83.9	
K4888-03	981-S-2- (20)	96	1.15	9.64	8.74	89.4	
K4889-01	DRUM-COMP	107	1.15	9.57	9.06	93.9	
K4893-01	RBR-200059-RT-1494	108	1.16	9.76	8.26	82.6	
K4897-01	RT-3425	109	1.12	9.73	8.92	90.6	
K4897-03	RBR-200028	110	1.14	9.64	9.14	94.1	
K4897-05	RBR-200052	111	1.13	9.65	9.37	96.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:27
Out Date: 09/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105112

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4898-01	FDS-1	112	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4898-02	WJC-1	113	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4898-03	6-7-EXJC-1	114	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4899-05	SVOC-GPC-BLANK	115	1.00	2.00	2.00	100.0	
K4899-06	PEST-GPC-BLANK	116	1.00	2.00	2.00	100.0	
K4899-07	PEST-GPC-BLANK-SPIKE	117	1.00	2.00	2.00	100.0	
K4899-08	PCB-GPC-BLANK	118	1.00	2.00	2.00	100.0	
K4899-09	PCB-GPC-BLANK-SPIKE	119	1.00	2.00	2.00	100.0	
K4899-10	SVOC-GPC2-BLANK	120	1.00	2.00	2.00	100.0	
K4899-11	PEST-GPC2-BLANK	121	1.00	2.00	2.00	100.0	
K4899-12	PEST-GPC2-BLANK-SPIKE	122	1.00	2.00	2.00	100.0	
K4899-13	PCB-GPC2-BLANK	123	1.00	2.00	2.00	100.0	
K4899-14	PCB-GPC2-BLANK -SPIKE	124	1.00	2.00	2.00	100.0	
K4900-01	RBR-200027-RBR-200049	125	1.14	9.94	9.58	95.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$