



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:10
In Date: 09/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 09/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105123

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4880-01	MH12T8	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-02	MH12T9	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-03	MH12W0	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-04	MH12W1	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-05	MH1305	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-06	MH1306	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-07	MH1307	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-08	MH1308	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-09	MH1309	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-10	MH1310	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-11	MH1311	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-12	MH1312	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-13	MH1313	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-14	MH1314	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-15	MH1315	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-16	MH1315D	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-17	MH1315S	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-18	MH1316	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-19	MH1317	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-20	MH1318	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K4880-21	MH1319	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$