



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/17/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:36
In Date: 09/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105165

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4666-12	PL-02-091619	1	1.15	9.71	9.37	96.0	
K4666-13	PL-02-091619-1	2	1.14	9.53	9.46	99.2	
K4666-14	PL-02-091619-2	3	1.18	9.69	9.27	95.1	
K4666-15	PL-02-091619-3	4	1.18	9.98	9.71	96.9	
K4666-16	PL-02-091619-4	5	1.13	9.7	9.39	96.4	
K4666-17	PL-02-091619-5	6	1.14	9.68	9.41	96.8	
K4910-01	390910009	7	1.13	9.55	1.68	6.5	SLUDGE SAMPLE
K4911-01	390905733	8	1.13	9.9	2.71	18.0	SLUDGE SAMPLE
K4912-01	390905734	9	1.14	9.97	2.92	20.2	SLUDGE SAMPLE
K4913-01	390905761	10	1.15	8.08	2.71	22.5	SLUDGE SAMPLE
K4914-01	390911704.2	11	1.14	9.98	1.31	1.9	SLUDGE SAMPLE
K4915-01	390912719	12	1.14	9.96	1.76	7.0	SLUDGE SAMPLE
K4916-01	390910091	13	1.14	9.97	5.88	53.7	SLUDGE SAMPLE
K4917-01	390910090	14	1.15	9.71	1.57	4.9	SLUDGE SAMPLE
K4918-01	SB-2	15	1.14	9.98	8.63	84.7	
K4918-02	SB-3	16	1.13	9.5	8.12	83.5	
K4918-03	SB-4	17	1.14	9.6	8.35	85.2	
K4918-04	SB-5	18	1.14	9.82	8.47	84.4	
K4918-05	SB-6	19	1.15	9.94	8.63	85.1	
K4918-06	SB-7	20	1.14	9.82	7.94	78.3	
K4918-07	SB-8	21	1.14	9.83	8.11	80.2	
K4918-08	SB-9	22	1.13	9.79	8.22	81.9	
K4918-09	SB-10	23	1.12	9.76	8.62	86.8	
K4918-10	SB-11	24	1.14	9.77	9.01	91.2	
K4918-11	SB-13	25	1.13	9.6	7.93	80.3	
K4918-12	SB-12	26	1.14	9.73	8.71	88.1	
K4918-13	SB-1	27	1.14	9.95	8.65	85.2	
K4920-01	1-2-WC-3	28	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/17/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:36
In Date: 09/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105165

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4920-02	6-7-EWC-3	29	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4920-03	FDS-1	30	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4920-04	DEXJC-1	31	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4920-06	WDC-1	32	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4920-07	WC-1	33	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100% SOLIDS
K4920-08	EXJC-1	34	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100 % SOLIDS
K4922-01	GRAB-40302	35	1.13	9.96	9.47	94.5	
K4923-01	SB-31-10.5-11.0	36	1.13	9.92	9.04	90.0	
K4923-02	SB-31-COMP	37	1.13	9.97	9.08	89.9	
K4923-03	SB-40-5.5-6.0	38	1.12	9.89	9.34	93.7	
K4923-04	SB-40-COMP	39	1.12	9.6	9.06	93.6	
K4923-05	SB-47-9.5-10.0	40	1.13	9.65	9.01	92.5	
K4923-06	SB-47-COMP	41	1.11	9.68	9.26	95.1	
K4923-07	SB-51-5.5-6.0	42	1.14	9.85	9.19	92.4	
K4923-08	SB-51-COMP	43	1.11	9.97	9.11	90.3	
K4923-09	SB-46-9.5-10.0	44	1.13	9.96	9.29	92.4	
K4923-10	SB-46-COMP	45	1.12	9.6	8.98	92.7	
K4924-01	SB-48-6.5-7.0	46	1.13	9.89	8.4	83.0	
K4924-02	SB-48-COMP	47	1.12	9.78	8.39	83.9	
K4924-03	SB-52-4.5-5.0	48	1.13	9.78	9.3	94.5	
K4924-04	SB-52-COMP	49	1.15	9.74	9.22	93.9	
K4924-05	SB-41-5.5-6.0	50	1.12	9.87	8.81	87.9	
K4924-06	SB-41-COMP	51	1.12	9.6	8.47	86.7	
K4924-07	SB-50-14.5-15.0	52	1.14	9.7	8.89	90.5	
K4924-08	SB-50-COMP	53	1.12	9.52	8.93	93.0	
K4924-09	SB-49-12.5-13.0	54	1.12	9.85	9.08	91.2	
K4924-10	SB-49-COMP	55	1.12	9.85	8.84	88.4	
K4924-11	SB-45-6.5-7.0	56	1.12	9.76	9.14	92.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/17/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:36
In Date: 09/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105165

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4924-12	SB-45-COMP	57	1.12	9.67	8.99	92.0	
K4924-13	SB-44-12.5-13.0	58	1.12	9.7	9.31	95.5	
K4924-14	SB-44-COMP	59	1.12	9.86	9.27	93.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$