



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:08
In Date: 09/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 09/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105265

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4665-01	EO-01-091919	59	1.13	9.68	9.05	92.6	
K4665-02	EO-01-091919	60	1.12	9.77	8.71	87.7	
K4665-03	EO-01-091919-1	61	1.12	9.51	8.39	86.7	
K4665-04	EO-01-091919-2	62	1.14	9.76	8.86	89.6	
K4665-05	EO-01-091919-3	63	1.13	9.93	9.09	90.5	
K4980-04	SAND-STOCKPILE-1	1	1.14	9.86	9.44	95.2	
K4980-05	SAND-STOCKPILE-1	2	1.16	9.82	9.4	95.2	
K4980-06	SAND-STOCKPILE-2	3	1.12	9.87	9.42	94.9	
K4980-07	SAND-STOCKPILE-2	4	1.16	9.97	9.28	92.2	
K4994-01	EXT-SW091-091819-2	5	1.16	9.94	9.3	92.7	
K4994-02	EXT-SW091-091819-3	6	1.14	9.84	9.07	91.1	
K4994-03	EXT-SW092-091819-1	7	1.14	9.82	9.22	93.1	
K4994-04	EXT-SW092-091819-2	8	1.13	9.86	8.5	84.4	
K4994-05	EXT-SW093-091819-1	9	1.14	9.76	8.74	88.2	
K4994-06	EXT-SW093-091819-2	10	1.16	9.61	8.72	89.5	
K4994-07	PST-B231-091919-1	11	1.14	9.85	9.18	92.3	
K4994-08	PST-B231-091919-2	12	1.13	9.53	9.00	93.7	
K4994-09	PST-SW117-091919-2	13	1.13	9.85	8.99	90.1	
K4994-10	PST-SW117-091919-3	14	1.14	9.86	9.03	90.5	
K4994-11	PST-B057-091919-1	15	1.13	9.67	8.84	90.3	
K4994-12	PST-B057-091919-2	16	1.13	9.65	8.5	86.5	
K4994-13	PST-B239-091919-1	17	1.12	9.94	8.99	89.2	
K4994-14	PST-B239-091919-2	18	1.12	9.61	8.7	89.3	
K4995-01	EXT-SW167-091819-1	19	1.12	9.54	8.54	88.1	
K4995-02	EXT-SW167-091819-2	20	1.12	9.62	8.76	89.9	
K4995-03	EXT-SW168-091819-1	21	1.12	9.5	8.57	88.9	
K4995-04	EXT-SW168-091819-2	22	1.13	9.6	8.6	88.2	
K4995-05	EXT-SW176-091819-1	23	1.13	9.7	8.97	91.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:08
In Date: 09/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 09/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105265

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4995-06	EXT-SW176-091819-2	24	1.12	9.78	9.1	92.1	
K4995-07	EXT-SW177-091819-1	25	1.15	9.77	9.22	93.6	
K4995-08	EXT-SW177-091819-2	26	1.14	9.96	9.41	93.8	
K4995-09	EXT-SW085-091819-3	27	1.13	9.65	9.15	94.1	
K4995-10	EXT-SW085-091819-4	28	1.15	9.55	9.08	94.4	
K4995-11	EXT-SW136-091819-1	29	1.15	9.66	9.43	97.3	
K4995-12	EXT-SW136-091819-2	30	1.13	9.63	9.03	92.9	
K4995-13	EXT-SW153-091819-1	31	1.14	9.63	9.06	93.3	
K4995-14	EXT-SW153-091819-2	32	1.14	9.75	9.22	93.8	
K4995-15	EXT-SW162-091819-1	33	1.14	9.77	9.42	95.9	
K4995-16	EXT-SW162-091819-2	34	1.15	9.71	9.24	94.5	
K4995-17	EXT-SW164-091819-1	35	1.14	9.66	9.26	95.3	
K4995-18	EXT-SW164-091819-2	36	1.14	9.77	9.28	94.3	
K4995-19	EXT-SW166-091819-1	37	1.15	9.65	8.76	89.5	
K4995-20	EXT-SW166-091819-2	38	1.12	9.87	9.36	94.2	
K4997-01	WC-08-091919	39	1.14	9.69	9.06	92.6	
K4997-02	WC-08-091919	40	1.14	9.78	9.15	92.7	
K4997-03	WC-09-091919	41	1.16	9.6	9.00	92.9	
K4997-04	WC-09-091919	42	1.16	9.67	9.06	92.8	
K4997-05	WC-10-091919	43	1.15	9.87	9.35	94.0	
K4997-06	WC-10-091919	44	1.12	9.61	9.1	94.0	
K4997-07	WC-11-091919	45	1.14	9.95	9.47	94.6	
K4997-08	WC-11-091919	46	1.14	9.73	9.31	95.1	
K4997-09	WC-12-091919	47	1.12	9.66	9.26	95.3	
K4997-10	WC-12-091919	48	1.13	9.63	9.16	94.5	
K5000-05	SVOC-GPC-BLANK	49	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-06	PEST-GPC-BLANK	50	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-07	PEST-GPC-BLANK-SPIKE	51	1.00	2.00	2.00	100.0	STORAGE BLK



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:08
In Date: 09/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 09/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105265

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5000-08	PCB-GPC-BLANK	52	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-09	PCB-GPC-BLANK-SPIKE	53	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-10	SVOC-GPC2-BLANK	54	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-11	PEST-GPC2-BLANK	55	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-12	PEST-GPC2-BLANK-SPIKE	56	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-13	PCB-GPC2-BLANK	57	1.00	2.00	2.00	100.0	STORAGE BLK
K5000-14	PCB-GPC2-BLANK -SPIKE	58	1.00	2.00	2.00	100.0	STORAGE BLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$