



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 09/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105332

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5046-01	MH1334	1	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-02	MH1335	2	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-03	MH1336	3	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-04	MH1337	4	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-05	MH1338	5	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-06	MH1338D	6	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-07	MH1338S	7	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-08	MH1339	8	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-09	MH1340	9	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-10	MH1341	10	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-11	MH1342	11	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-12	MH1343	12	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-13	MH1344	13	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-14	MH1345	14	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-15	MH1346	15	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-16	MH1348	16	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-17	MH1349	17	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-18	MH1350	18	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-19	MH1351	19	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-20	MH1352	20	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-21	MH1353	21	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS
K5046-22	MH1354	22	1.00	2.00	2.00	100.0	DUST SAMPLE 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$