



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:10
In Date: 09/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105335

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5050-01	JLMB4	1	1.14	9.82	8.14	80.6	
K5050-02	JLMB5	2	1.16	9.79	6.9	66.5	
K5050-03	JLMB6	3	1.16	9.82	8.12	80.4	
K5050-04	JLMB7	4	1.15	9.94	7.7	74.5	
K5050-05	JLMB9	5	1.12	9.6	8.31	84.8	
K5050-06	JLMC0	6	1.12	9.73	8.07	80.7	
K5050-07	JLMC1	7	1.18	9.89	7.71	75.0	
K5050-08	JLMC4	8	1.14	9.84	7.92	77.9	
K5050-09	JLMC5	9	1.19	9.98	6.75	63.3	
K5050-10	JLMC6	10	1.16	9.6	7.51	75.2	
K5050-11	JLMC9	11	1.12	9.54	7.36	74.1	
K5050-12	JLMD0	12	1.14	9.94	7.37	70.8	
K5050-13	JLMD1	13	1.13	9.54	7.74	78.6	
K5050-14	JLMD2	14	1.17	9.92	7.69	74.5	
K5050-15	JLMD3	15	1.15	9.98	7.68	74.0	
K5050-16	JLMD4	16	1.1	9.93	7.56	73.2	
K5050-17	JLMD5	17	1.17	9.96	7.86	76.1	
K5050-18	JLMD6	18	1.18	9.97	7.48	71.7	
K5050-19	JLMD7	19	1.16	9.96	7.9	76.6	
K5050-20	JLME0	20	1.13	9.96	8.6	84.6	
K5050-21	JLME0MS	21	1.13	9.96	8.6	84.6	
K5050-22	JLME0MSD	22	1.13	9.96	8.6	84.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$