



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:10
In Date: 09/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105336

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5051-01	MJLMB4	1	1.14	9.82	8.14	80.6	
K5051-02	MJLMB5	2	1.16	9.79	6.9	66.5	
K5051-03	MJLMB6	3	1.16	9.82	8.12	80.4	
K5051-04	MJLMB7	4	1.15	9.94	7.7	74.5	
K5051-05	MJLMB9	5	1.12	9.6	8.31	84.8	
K5051-06	MJLMC0	6	1.12	9.73	8.07	80.7	
K5051-07	MJLMC1	7	1.18	9.89	7.71	75.0	
K5051-08	MJLMC4	8	1.14	9.84	7.92	77.9	
K5051-09	MJLMC5	9	1.19	9.98	6.75	63.3	
K5051-10	MJLMC6	10	1.16	9.6	7.51	75.2	
K5051-11	MJLMC9	11	1.12	9.54	7.36	74.1	
K5051-12	MJLMD0	12	1.14	9.94	7.37	70.8	
K5051-13	MJLMD1	13	1.13	9.54	7.74	78.6	
K5051-14	MJLMD2	14	1.17	9.92	7.69	74.5	
K5051-15	MJLMD3	15	1.15	9.98	7.68	74.0	
K5051-16	MJLMD4	16	1.1	9.93	7.56	73.2	
K5051-17	MJLMD5	17	1.17	9.96	7.86	76.1	
K5051-18	MJLMD6	18	1.18	9.97	7.48	71.7	
K5051-19	MJLMD7	19	1.16	9.96	7.9	76.6	
K5051-20	MJLME0	20	1.13	9.96	8.6	84.6	
K5051-21	MJLME0D	21	1.13	9.96	8.6	84.6	
K5051-22	MJLME0S	22	1.13	9.96	8.6	84.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$