



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/27/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105340

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K4668-16	SU-04-092619	12	1.15	9.58	8.98	92.9	
K4668-17	SU-04-092619-1	13	1.15	9.85	8.41	83.4	
K4668-18	SU-04-092619-2	14	1.13	9.6	8.76	90.1	
K4822-03	SB-3A- (15-16)	1	1.15	9.91	8.3	81.6	
K4822-12	SB-3B- (15-16)	2	1.14	9.61	7.92	80.0	
K5041-01	TH-2	3	1.13	9.81	9.00	90.7	
K5041-02	TH-3	4	1.13	9.6	8.66	88.9	
K5041-03	TH-4	5	1.14	9.77	8.93	90.3	
K5041-04	TH-5	6	1.12	9.73	8.51	85.8	
K5041-05	TP-1-5	7	1.12	9.91	8.9	88.5	
K5055-01	1-STAIRWELL	8	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100 % SOLIDS
K5055-02	2-STAIRWELL	9	1.00	2.00	2.00	100.0	CAULKING SAMPLE 100 % SOLIDS
K5060-01	FK-PS-01-060	10	1.15	9.97	8.34	81.5	
K5060-02	FK-PS-01-060	11	1.12	9.69	8.2	82.6	
K5064-01	DRUM-COMP	15	1.14	9.72	9.47	97.1	
K5067-01	ROLL-OFF-1	16	1.14	9.73	8.55	86.3	
K5067-02	ROLL-OFF-2	17	1.13	9.89	8.14	80.0	
K5067-03	ROLL-OFF-3	18	1.14	9.68	7.93	79.5	
K5067-04	ROLL-OFF-4	19	1.14	9.54	8.31	85.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$