



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105350

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5072-01	0689-SED01	1	1.14	9.9	8.09	79.3	
K5072-02	K5072-01MS	2	1.14	9.9	8.09	79.3	
K5072-04	0689-SED02	3	1.16	9.83	8.02	79.1	
K5072-05	0689-SED03	4	1.15	9.6	4.15	35.5	
K5072-06	0689-SED04	5	1.14	9.65	6.48	62.7	
K5072-07	0689-SED05	6	1.15	9.85	2.48	15.3	
K5072-08	0689-SED06	7	1.15	9.6	2.53	16.3	
K5072-09	0689-SED07	8	1.15	9.86	2.78	18.7	
K5072-10	0689-SED08	9	1.14	9.93	2.63	17.0	
K5072-11	0689-SED09	10	1.17	9.57	4.44	38.9	
K5072-12	0689-SED10	11	1.15	9.82	6.25	58.8	
K5073-01	P002-SD001-0006-01	12	1.15	9.51	6.17	60.0	
K5073-02	P002-SD001-1824-01	13	1.15	9.71	3.06	22.3	
K5073-03	P002-SD002-0006-01	14	1.15	9.93	7.68	74.4	
K5073-04	P002-SD002-1824-01	15	1.14	9.61	2.99	21.8	
K5073-05	P002-SD003-0006-01	16	1.14	9.97	8.53	83.7	
K5073-06	P002-SD003-1824-01	17	1.14	9.62	2.9	20.8	
K5073-07	P002-SD004-0006-01	18	1.14	9.52	4.07	35.0	
K5073-08	K5073-07MS	19	1.14	9.52	4.07	35.0	
K5073-10	P002-SD004-0006-02	20	1.14	9.52	2.99	22.1	
K5073-11	P002-SD004-1824-01	21	1.15	9.73	7.1	69.3	
K5073-12	P002-SD005-0006-01	22	1.14	9.77	3.23	24.2	
K5073-13	P002-SD005-1824-01	23	1.14	9.5	2.74	19.1	
K5073-14	P002-SD006-0006-01	24	1.14	9.61	4.28	37.1	
K5073-15	P002-SD006-1824-01	25	1.17	9.83	8.16	80.7	
K5073-16	P002-SD007-0006-01	26	1.15	9.55	2.91	21.0	
K5073-17	P002-SD007-1824-01	27	1.16	9.8	6.2	58.3	
K5073-18	P002-SD008-0006-01	28	1.15	9.64	3.18	23.9	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105350

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5073-19	P002-SD008-1824-01	29	1.15	9.95	4.68	40.1	
K5073-20	P002-SD009-0006-01	30	1.15	9.63	3.18	23.9	
K5073-21	P002-SD009-1824-01	31	1.15	9.58	2.78	19.3	
K5073-22	P002-SD010-0006-01	32	1.13	9.73	3.13	23.3	
K5074-01	P002-SD010-1824-01	33	1.14	9.57	2.76	19.2	
K5074-02	P002-SD011-0006-01	34	1.14	9.58	4.59	40.9	
K5074-03	P002-SD011-1824-01	35	1.15	9.79	6.43	61.1	
K5074-04	P002-SD012-0006-01	36	1.15	9.76	2.51	15.8	
K5074-05	P002-SD012-1824-01	37	1.15	9.91	2.19	11.9	
K5074-06	P003-SD001-0006-01	38	1.14	9.71	2.7	18.2	
K5074-07	P003-SD001-1824-01	39	1.15	9.56	2.18	12.2	
K5074-08	P003-SD002-0006-01	40	1.16	9.53	4.63	41.5	
K5074-09	P003-SD002-1824-01	41	1.15	9.87	8.15	80.3	
K5074-10	K5074-09MS	42	1.15	9.87	8.15	80.3	
K5074-12	P003-SD002-1824-02	43	1.15	9.95	8.04	78.3	
K5074-13	P003-SD003-0006-01	44	1.15	9.84	2.17	11.7	
K5074-14	P003-SD003-1824-01	45	1.15	9.74	2.00	9.9	
K5075-01	0689-SED11	46	1.15	9.66	7.4	73.4	
K5075-02	0689-SED12	47	1.17	9.93	7.8	75.7	
K5075-03	P001-SD001-0006-01	48	1.14	9.64	9.25	95.4	
K5075-04	P001-SD001-1824-01	49	1.13	9.74	4.16	35.2	
K5075-05	P001-SD002-0006-01	50	1.15	9.56	7.71	78.0	
K5075-06	P001-SD002-1824-01	51	1.13	9.66	8.39	85.1	
K5075-07	P001-SD003-0006-01	52	1.14	9.87	4.48	38.3	
K5075-08	P001-SD003-1824-01	53	1.14	9.65	7.86	79.0	
K5075-09	P001-SD004-0006-01	54	1.14	9.7	3.73	30.3	
K5075-10	P001-SD004-1824-01	55	1.13	9.84	2.64	17.3	
K5075-11	P001-SD005-0006-01	56	1.14	9.8	3.02	21.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105350

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5075-12	P001-SD005-1824-01	57	1.14	9.94	4.74	40.9	
K5075-13	P001-SD006-0006-01	58	1.13	9.53	3.00	22.3	
K5075-14	P001-SD006-1824-01	59	1.14	9.65	7.04	69.3	
K5075-15	K5075-14MS	60	1.14	9.65	7.04	69.3	
K5075-17	P001-SD006-1824-02	61	1.13	9.81	7.43	72.6	
K5077-01	PST-024-B271-092519	62	1.15	9.81	8.65	86.6	
K5077-02	PST-024-B272-092519	63	1.14	9.5	8.88	92.6	
K5077-03	PST-024-B273-092519	64	1.14	9.74	8.57	86.4	
K5077-04	PST-024-B274-092519	65	1.16	9.55	9.26	96.5	
K5077-05	PST-024-B275-092519	66	1.14	9.8	9.08	91.7	
K5077-06	PST-024-B276-092519	67	1.15	9.68	9.28	95.3	
K5077-07	PST-024-B277-092519	68	1.15	9.71	8.79	89.3	
K5077-08	PST-024-B278-092519	69	1.13	9.86	9.61	97.1	
K5077-09	PST-024-B279-092519	70	1.15	9.52	9.07	94.6	
K5077-10	PST-024-B280-092519	71	1.12	9.75	9.03	91.7	
K5077-11	EXT-024-SW075-092519	72	1.15	9.7	9.16	93.7	
K5077-12	PST-024-B262-092519	73	1.16	9.72	9.37	95.9	
K5077-13	PST-024-B263-092519	74	1.13	9.83	7.96	78.5	
K5077-14	PST-024-B264-092619	75	1.13	9.58	8.1	82.5	
K5077-15	PST-024-B265-092519	76	1.14	9.53	9.15	95.5	
K5077-16	PST-024-B266-092519	77	1.15	9.54	8.96	93.1	
K5077-17	PST-024-B267-092519	78	1.14	9.64	9.3	96.0	
K5077-18	PST-024-B268-092519	79	1.13	9.63	9.08	93.5	
K5077-19	PST-024-B269-092519	80	1.13	9.71	9.06	92.4	
K5077-20	PST-024-B270-092619	81	1.12	9.85	9.02	90.5	
K5077-21	PST-024-SW136-092619	82	1.13	9.7	9.44	97.0	
K5077-22	PST-024-SW137-092519	83	1.16	9.6	9.22	95.5	
K5077-23	PST-024-SW138-092519	84	1.14	9.87	9.35	94.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105350

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5078-01	EXT-024-SW065-092519	85	1.14	9.63	9.06	93.3	
K5078-02	EXT-024-SW066-092519	86	1.13	9.56	8.5	87.4	
K5078-03	EXT-024-SW067-092519	87	1.13	9.9	9.66	97.3	
K5078-04	EXT-024-SW068-092619	88	1.16	9.56	8.81	91.1	
K5078-05	EXT-024-SW069-092519	89	1.15	9.77	9.41	95.8	
K5078-06	EXT-024-SW070-092519	90	1.15	9.63	9.18	94.7	
K5078-07	EXT-024-SW071-092619	91	1.16	9.96	9.61	96.0	
K5078-08	EXT-024-SW072-092619	92	1.16	9.92	9.62	96.6	
K5078-09	EXT-024-SW073-092519	93	1.15	9.65	9.29	95.8	
K5078-10	EXT-024-SW074-092519	94	1.15	9.54	8.92	92.6	
K5078-11	PST-024-B281-092519	95	1.15	9.9	9.56	96.1	
K5078-12	PST-024-B282-092619	96	1.16	9.61	9.22	95.4	
K5078-13	PST-024-SW128-092519	97	1.15	9.63	9.1	93.7	
K5078-14	PST-024-SW129-092519	98	1.15	9.6	8.85	91.1	
K5078-15	PST-024-SW130-092519	99	1.16	9.55	8.14	83.2	
K5078-16	PST-024-SW131-092519	100	1.15	9.59	8.17	83.2	
K5078-17	PST-024-SW132-092519	101	1.13	9.92	8.98	89.3	
K5078-18	PST-024-SW133-092519	102	1.12	9.66	9.26	95.3	
K5078-19	PST-024-SW134-092519	103	1.13	9.81	8.62	86.3	
K5078-20	PST-024-SW135-092519	104	1.15	9.63	9.12	94.0	
K5083-01	20190708-RR-COMPOSITE-8	105	1.14	9.86	5.98	55.5	
K5083-02	20190707-RR-COMPOSITE-7	106	1.14	9.74	6.01	56.6	
K5083-03	20190706-RR-COMP-6	107	1.16	9.58	6.11	58.8	
K5083-04	20190705-RR-COMPOSITE-5	108	1.13	9.67	6.33	60.9	
K5083-05	20190704-RR-COMPOSITE-4	109	1.16	9.52	7.04	70.3	
K5083-06	20190703-RR-COMPOSITE-3	110	1.15	9.54	6.53	64.1	
K5083-07	20190702-RR-COMPOSITE-2	111	1.15	9.77	6.11	57.5	
K5083-08	20190701-RR-COMPOSITE-1	112	1.13	9.86	5.88	54.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105350

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5087-05	SVOC-GPC-BLANK	113	1.00	2.00	2.00	100.0	
K5087-06	PEST-GPC-BLANK	114	1.00	2.00	2.00	100.0	
K5087-07	PEST-GPC-BLANK-SPIKE	115	1.00	2.00	2.00	100.0	
K5087-08	PCB-GPC-BLANK	116	1.00	2.00	2.00	100.0	
K5087-09	PCB-GPC-BLANK-SPIKE	117	1.00	2.00	2.00	100.0	
K5087-10	SVOC-GPC2-BLANK	118	1.00	2.00	2.00	100.0	
K5087-11	PEST-GPC2-BLANK	119	1.00	2.00	2.00	100.0	
K5087-12	PEST-GPC2-BLANK-SPIKE	120	1.00	2.00	2.00	100.0	
K5087-13	PCB-GPC2-BLANK	121	1.00	2.00	2.00	100.0	
K5087-14	PCB-GPC2-BLANK -SPIKE	122	1.00	2.00	2.00	100.0	
K5093-01	SB-01-(9.5-10)	123	1.14	9.97	9.32	92.6	
K5093-02	SB-01-COMP	124	1.17	9.71	9.31	95.3	
K5093-03	SB-02-(9.5-10)	125	1.16	9.77	8.53	85.6	
K5093-04	SB-02-COMP	126	1.12	9.94	8.71	86.1	
K5093-05	SB-03-(9.5-10)	127	1.14	9.72	8.62	87.2	
K5093-06	SB-03-COMP	128	1.14	9.52	8.82	91.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$