



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 13:09
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105352

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5026-01	C0AW4	1	1.16	9.7	7.23	71.1	
K5026-02	C0AW5	2	1.14	9.87	7.26	70.1	
K5026-03	C0AW6	3	1.14	9.73	6.71	64.8	
K5026-04	C0AW7	4	1.14	9.73	8.66	87.5	
K5026-05	C0AW8	5	1.12	9.66	8.45	85.8	
K5026-06	C0AW9	6	1.14	9.52	8.89	92.5	
K5026-09	C0AX4	7	1.12	9.85	7.56	73.8	
K5026-10	C0AX5	8	1.12	9.81	7.73	76.1	
K5026-11	C0AX6	9	1.14	9.68	8.34	84.3	
K5026-12	C0AX7	10	1.14	9.63	8.00	80.8	
K5026-13	C0AX8	11	1.12	9.85	8.74	87.3	
K5026-14	C0B02	12	1.13	9.65	7.08	69.8	
K5026-15	C0B03	13	1.13	9.72	7.42	73.2	
K5026-16	C0B04	14	1.13	9.68	7.58	75.4	
K5026-17	C0B05	15	1.13	9.7	8.42	85.1	
K5026-18	C0B06	16	1.13	9.62	8.85	90.9	
K5026-19	C0B11	17	1.14	9.5	8.62	89.5	
K5026-20	C0B11MS	18	1.14	9.5	8.62	89.5	
K5026-21	C0B11MSD	19	1.14	9.5	8.62	89.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$