



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:30
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105360

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5056-01	DBAC7	1	1.00	2.00	2.00	100.0	P.E.SAMPLE
K5056-03	DBAC9	2	1.00	2.00	2.00	100.0	P.E.SAMPLE
K5056-04	DBAD8	3	1.15	9.88	6.71	63.7	
K5056-05	DBAE0	4	1.12	9.76	8.48	85.2	
K5056-06	DBAE1	5	1.16	9.81	7.22	70.1	
K5056-07	DBAE2	6	1.18	9.86	7.06	67.7	
K5056-08	DBAE3	7	1.15	9.94	8.2	80.2	
K5056-09	DBAE5	8	1.15	9.96	7.17	68.3	
K5056-10	DBAE6	9	1.15	9.71	7.94	79.3	
K5056-11	DBAE9	10	1.11	9.72	8.29	83.4	
K5056-12	DBAE9MS	11	1.11	9.72	8.29	83.4	
K5056-13	DBAE9MSD	12	1.11	9.72	8.29	83.4	
K5056-14	DBAF0	13	1.16	9.58	8.3	84.8	
K5056-15	DBAF1	14	1.14	9.87	9.00	90.0	
K5056-16	DBAF2	15	1.16	9.77	7.93	78.6	
K5056-17	DBAF3	16	1.15	9.75	8.43	84.7	
K5056-18	DBAF4	17	1.14	9.62	8.57	87.6	
K5056-19	DBAF5	18	1.13	9.54	8.85	91.8	
K5056-20	DBAF6	19	1.19	9.82	8.19	81.1	
K5056-21	DBAG5	20	1.13	9.65	6.68	65.1	
K5056-22	DBAG6	21	1.16	9.97	7.4	70.8	
K5056-23	DBAG7	22	1.16	9.94	8.65	85.3	
K5056-24	DBAH0	23	1.14	9.74	6.11	57.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$