



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 9/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:05
In Date: 09/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105362

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5063-01	MDBAD6	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
K5063-02	MDBAH8	2	1.16	9.82	8.03	79.3	
K5063-03	MDBAJ0	3	1.13	9.87	7.91	77.6	
K5063-04	MDBAJ2	4	1.14	9.66	8.74	89.2	
K5063-05	MDBAJ3	5	1.12	9.8	7.59	74.5	
K5063-06	MDBAJ5	6	1.14	9.97	8.08	78.6	
K5063-07	MDBAJ6	7	1.15	9.84	7.86	77.2	
K5063-08	MDBAJ9	8	1.14	9.75	7.22	70.6	
K5063-09	MDBAK2	9	1.15	9.83	3.62	28.5	
K5063-10	MDBAK2D	10	1.15	9.83	3.62	28.5	
K5063-11	MDBAK2S	11	1.15	9.83	3.62	28.5	
K5063-12	MDBAK3	12	1.16	9.8	3.37	25.6	
K5063-13	MDBAK4	13	1.14	9.75	3.98	33.0	
K5063-14	MDBAK5	14	1.14	9.96	2.63	16.9	
K5063-15	MDBAK6	15	1.14	9.98	2.59	16.4	
K5063-16	MDBAK7	16	1.15	9.97	5.31	47.2	
K5063-17	MDBAK9	17	1.16	9.67	5.27	48.3	
K5063-18	MDBAL1	18	1.14	9.87	6.56	62.1	
K5063-19	MDBAF3	19	1.17	9.66	8.29	83.9	
K5063-20	MDBAF4	20	1.16	9.7	8.64	87.6	
K5063-21	MDBAF5	21	1.14	9.97	8.26	80.6	
K5063-22	MDBAK8	22	1.16	9.67	8.93	91.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$