



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/1/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 09/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 10/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105382

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5066-01	MDBAF7	1	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-02	MDBAF8	2	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-03	MDBAF9	3	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-04	MDBAF9D	4	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-05	MDBAF9S	5	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-06	MDBAG0	6	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-07	MDBAG2	7	1.00	2.00	2.00	100.0	100 % SOLIDS
K5066-08	MDBAG4	8	1.00	2.00	2.00	100.0	100 % SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$