



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 13:47
In Date: 10/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 10/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105403

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5098-01	MJLMA2	1	1.14	9.56	4.63	41.4	
K5098-02	MJLMA3	2	1.15	9.98	5.04	44.1	
K5098-03	MJLMA4	3	1.16	9.78	7.03	68.1	
K5098-04	MJLMA5	4	1.13	9.65	6.54	63.5	
K5098-05	MJLMD9	5	1.12	9.6	8.12	82.5	
K5098-06	MJLME1	6	1.12	9.55	7.82	79.5	
K5098-07	MJLME2	7	1.13	9.82	7.57	74.1	
K5098-08	MJLME4	8	1.14	9.56	7.64	77.2	
K5098-09	MJLME5	9	1.14	9.54	7.14	71.4	
K5098-10	MJLME6	10	1.13	9.75	7.71	76.3	
K5098-11	MJLME7	11	1.14	9.58	7.23	72.2	
K5098-12	MJLME9	12	1.12	9.61	7.26	72.3	
K5098-13	MJLMF0	13	1.12	9.52	7.93	81.1	
K5098-14	MJLMF1	14	1.14	9.63	7.12	70.4	
K5098-15	MJLMF2	15	1.12	9.58	7.35	73.6	
K5098-16	MJLMF4	16	1.12	9.53	7.56	76.6	
K5098-17	MJLMF5	17	1.11	9.8	7.39	72.3	
K5098-18	MJLMF6	18	1.13	9.98	7.88	76.3	
K5098-19	MJLMF7	19	1.12	9.55	7.31	73.4	
K5098-20	MJLM82	20	1.12	9.72	7.96	79.5	
K5098-21	MJLM82D	21	1.12	9.72	7.96	79.5	
K5098-22	MJLM82S	22	1.12	9.72	7.96	79.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$