



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: JIGNESH

Date: 10/7/2019

OVENTEMP IN Celsius(°C): 107

Time IN: 15:13

In Date: 10/04/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 07:35

Out Date: 10/05/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB105505

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
K5126-01	MJLM50	1	1.13	9.82	8.57	85.6	
K5126-02	MJLM51	2	1.14	9.7	8.64	87.6	
K5126-03	MJLM54	3	1.13	9.98	6.03	55.4	
K5126-04	MJLM66	4	1.14	9.57	7.3	73.1	
K5126-05	MJLM67	5	1.14	9.77	6.78	65.4	
K5126-06	MJLM74	6	1.14	9.75	6.96	67.6	
K5126-07	MJLM75	7	1.15	9.67	6.71	65.3	
K5126-08	MJLM78	8	1.16	9.83	4.58	39.4	
K5126-09	MJLM79	9	1.15	9.6	7.95	80.5	
K5126-10	MJLM80	10	1.14	9.95	6.68	62.9	
K5126-11	MJLM83	11	1.15	9.98	7.26	69.2	
K5126-12	MJLM83D	12	1.15	9.98	7.26	69.2	
K5126-13	MJLM83S	13	1.15	9.98	7.26	69.2	
K5126-14	MJLM90	14	1.14	9.96	8.97	88.8	
K5126-15	MJLM91	15	1.16	9.83	8.51	84.8	
K5126-16	MJLM92	16	1.16	9.97	7.36	70.4	
K5126-17	MJLMA6	17	1.14	9.68	4.89	43.9	
K5126-18	MJLMA7	18	1.14	9.98	6.45	60.1	
K5126-19	MJLMA8	19	1.17	9.77	4.55	39.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$