



PERCENT SOLID

Supervisor: AMANDEEP
Analyst: JIGNESH
Date: 10/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:13
In Date: 10/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 10/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105514

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5124-01	JLM50	1	1.13	9.82	8.57	85.6	
K5124-02	JLM51	2	1.14	9.7	8.64	87.6	
K5124-03	JLM54	3	1.13	9.98	6.03	55.4	
K5124-04	JLM66	4	1.14	9.57	7.3	73.1	
K5124-05	JLM67	5	1.14	9.77	6.78	65.4	
K5124-06	JLM74	6	1.14	9.75	6.96	67.6	
K5124-07	JLM75	7	1.15	9.67	6.71	65.3	
K5124-08	JLM78	8	1.16	9.83	4.58	39.4	
K5124-09	JLM79	9	1.15	9.6	7.95	80.5	
K5124-10	JLM80	10	1.14	9.95	6.68	62.9	
K5124-11	JLM83	11	1.15	9.98	7.26	69.2	
K5124-12	JLM83MS	12	1.15	9.98	7.26	69.2	
K5124-13	JLM83MSD	13	1.15	9.98	7.26	69.2	
K5124-14	JLM90	14	1.14	9.96	8.97	88.8	
K5124-15	JLM91	15	1.16	9.83	8.51	84.8	
K5124-16	JLM92	16	1.16	9.97	7.36	70.4	
K5124-17	JLMA6	17	1.14	9.68	4.89	43.9	
K5124-18	JLMA7	18	1.14	9.98	6.45	60.1	
K5124-19	JLMA8	19	1.17	9.77	4.55	39.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$