



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:17
In Date: 10/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:49
Out Date: 10/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105576

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5272-01	MGAPP7	1	1.14	9.6	7.65	77.0	
K5272-02	MGAPP8	2	1.12	9.98	7.96	77.2	
K5272-03	MGAPP9	3	1.12	9.9	9.57	96.2	
K5272-04	MGAPP9D	4	1.12	9.9	9.57	96.2	
K5272-05	MGAPP9S	5	1.12	9.9	9.57	96.2	
K5272-06	MGAPQ0	6	1.13	9.78	8.78	88.4	
K5272-07	MGAPQ1	7	1.13	9.62	8.92	91.8	
K5272-08	MGAPQ2	8	1.13	9.63	9.04	93.1	
K5272-09	MGAPQ3	9	1.13	9.67	7.81	78.2	
K5272-10	MGAPQ4	10	1.13	9.59	7.15	71.2	
K5272-11	MGAPQ5	11	1.13	9.97	7.48	71.8	
K5272-12	MGAPQ6	12	1.14	9.72	8.46	85.3	
K5272-13	MGAPQ7	13	1.12	9.8	8.53	85.4	
K5272-14	MGAPQ8	14	1.14	9.72	8.73	88.5	
K5272-15	MGAPQ9	15	1.15	9.94	8.64	85.2	
K5272-16	MGAPR0	16	1.15	9.98	8.81	86.7	
K5272-17	MGAPR1	17	1.15	9.65	9.33	96.2	
K5272-18	MGAPR2	18	1.13	9.8	9.39	95.3	
K5272-19	MGAPR3	19	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$