



PERCENT SOLID

Supervisor: apatel  
Analyst: JIGNESH  
Date: 10/11/2019

OVENTEMP IN Celsius(°C): 107  
Time IN: 14:32  
In Date: 10/10/2019  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 07:25  
Out Date: 10/11/2019  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: %SOLIDS-OVEN

QC:LB105614

| Lab ID   | Client SampleID | Dish # | Dish Wt (g) (A) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments   |
|----------|-----------------|--------|-----------------|--------------------------|----------------------------|---------|------------|
| K5308-06 | X2J24           | 1      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E.SAMPLE |
| K5308-07 | X2J25           | 2      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E.SAMPLE |
| K5308-08 | X2J26           | 3      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E.SAMPLE |
| K5308-09 | X2J27           | 4      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E.SAMPLE |
| K5308-10 | X2J28           | 5      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E.SAMPLE |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$