



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 10/14/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 10/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:27
Out Date: 10/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB105640

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K5345-01	ETOJ6	1	1.13	9.58	8.27	84.5	
K5345-02	ETOJ6MS	2	1.13	9.58	8.27	84.5	
K5345-03	ETOJ6MSD	3	1.13	9.58	8.27	84.5	
K5345-04	ETOJ7	4	1.14	9.53	7.98	81.5	
K5345-05	ETOJ8	5	1.13	9.75	8.03	80.0	
K5345-06	ETOJ9	6	1.12	9.64	8.23	83.5	
K5345-07	ETOK0	7	1.12	9.7	8.6	87.2	
K5345-08	ETOK1	8	1.13	9.66	8.48	86.2	
K5345-09	ETOK2	9	1.12	9.66	8.36	84.8	
K5345-10	ETOK3	10	1.14	9.67	8.24	83.2	
K5345-11	ETOK4	11	1.14	9.8	8.19	81.4	
K5345-12	ETOK5	12	1.14	9.85	8.11	80.0	
K5345-13	ETOK6	13	1.15	9.88	8.43	83.4	
K5345-14	ETOK7	14	1.14	9.91	8.7	86.2	
K5345-15	ETOK8	15	1.15	9.56	8.34	85.5	
K5345-16	ETOK9	16	1.12	9.6	8.16	83.0	
K5345-17	ETOL0	17	1.13	9.73	8.02	80.1	
K5345-18	ETOL2	18	1.15	9.65	8.4	85.3	
K5345-19	ETOL3	19	1.12	9.88	8.61	85.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$